

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | pb97604       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 08 Min 06 Sec   |
| Total Counter :  | 600250        | Total Running Time : | 2275 Hour 44 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 48                | %       | Manual ;           | 57           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6400000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 35                | Min     | Supply Time ;      | 1.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 9.74              | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 9.35              | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 7.95              | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 7.79              | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 7.64              | kN      | Yes 15             | Times        |           |
| LEP :                              | 20.0 %      | 6.23              | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 5.85              | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |            |      |       |            |    |
|----------------------|---------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 1.58          | kN    | Min :      | 1.44 | kN    | Punch NO.; | 26 |
| SD :                 | 3.10          | %     | Max :      | 1.67 | kN    | Punch NO.; | 37 |
| Mean Main Pressure : | 7.87          | kN    | Min :      | 6.55 | kN    | Punch NO.; | 26 |
| SD :                 | 7.30          | %     | Max :      | 8.84 | kN    | Punch NO.; | 10 |
| Filling Depth :      | Increase ; 26 | Times | Decrease ; | 66   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.37    | mm  | Main Thickness ; | 1.37              | mm           |
| Filling Depth :               | 3.57    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.818   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.32 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 34 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 600250            | Tablets |       |             |         |
| A.W.C Counter : | 587003            | Tablets |       |             |         |
| HEP :           | 375 (0.1%)        | Tablets | LEP : | 1789 (0.3%) | Tablets |
| GOOD :          | 584839<br>(99.6%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | pb97604       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 08 Min 06 Sec   |
| Total Counter :  | 600250        | Total Running Time : | 2275 Hour 44 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 46                | %       | Manual ;           | 30    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 35                | Min     | Supply Time ;      | 1.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 9.84              | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 9.44              | kN      | Yes                | 15    | Rotations 20 Tablets |
| HCP :                              | 2.0 %  | 8.03              | kN      | Yes                | 20    | Times                |
| Ref :                              |        | 7.87              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 7.71              | kN      | Yes                | 20    | Times                |
| LEP :                              | 20.0 % | 6.29              | kN      | Yes                | 15    | Rotations 20 Tablets |
| LSP :                              | 25.0 % | 5.90              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.80       | kN | Min : | 1.73       | kN | Punch NO.; | 12 |
| SD :                 | 2.30       | %  | Max : | 1.91       | kN | Punch NO.; | 2  |
| Mean Main Pressure : | 7.88       | kN | Min : | 6.76       | kN | Punch NO.; | 55 |
| SD :                 | 5.70       | %  | Max : | 9.30       | kN | Punch NO.; | 10 |
| Filling Depth :      | Increase ; | 14 | Times | Decrease ; | 44 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.57    | mm  | Main Thickness ; | 1.57              | mm           |
| Filling Depth :               | 3.99    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.818   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.33 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 34 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 600250            | Tablets |       |             |         |
| A.W.C Counter : | 587003            | Tablets |       |             |         |
| HEP :           | 375 (0.1%)        | Tablets | LEP : | 1789 (0.3%) | Tablets |
| GOOD :          | 584839<br>(99.6%) | Tablets |       |             |         |