

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | sj97904       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 06 Min 00 Sec   |
| Total Counter :  | 6251218       | Total Running Time : | 2304 Hour 44 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 36                | %       | Manual ;           | 57    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 90                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 190               | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 9.73              | kN      | No                 |       |                      |
| HEP :                              | 22.0 % | 9.13              | kN      | Yes                | 10    | Rotations 20 Tablets |
| HCP :                              | 2.5 %  | 7.67              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 7.49              | kN      |                    |       |                      |
| LCP :                              | 2.5 %  | 7.30              | kN      | Yes                | 15    | Times                |
| LEP :                              | 28.0 % | 5.39              | kN      | Yes                | 15    | Rotations 20 Tablets |
| LSP :                              | 35.0 % | 4.87              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |      |       |            |    |
|----------------------|------------|-------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.47  | kN    | Min :      | 1.23 | kN    | Punch NO.; | 19 |
|                      |            | 8.60  | %     | Max :      | 1.72 | kN    | Punch NO.; | 37 |
| Mean Main Pressure : | SD :       | 4.19  | kN    | Min :      | 1.94 | kN    | Punch NO.; | 61 |
|                      |            | 35.60 | %     | Max :      | 7.01 | kN    | Punch NO.; | 13 |
| Filling Depth :      | Increase ; | 163   | Times | Decrease ; | 225  | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 0        | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.40    | mm  | Pre Thickness ;  | 1.40     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.37    | mm  | Main Thickness ; | 1.37     | mm           |
| Filling Depth :               | 4.28    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.821   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 68      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 152     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | On      |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.34 | kN |             |    |
|                  | Max ;     | 0.38 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                 |         |       |              |         |
|-----------------|-----------------|---------|-------|--------------|---------|
| Total Counter : | 6251218         | Tablets |       |              |         |
| A.W.C Counter : | 6217796         | Tablets |       |              |         |
| HEP :           | 2436 (0.0%)     | Tablets | LEP : | 64897 (1.0%) | Tablets |
| GOOD :          | 6150463 (98.9%) | Tablets |       |              |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PROPRA 10MG<br>STPJ2A00 |
| User ID :        | sj97904       | Batch NO. :          | PJM2500610              |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 06 Min 00 Sec   |
| Total Counter :  | 6251218       | Total Running Time : | 2304 Hour 44 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 38                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6400000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 90                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 190               | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 10.47             | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 10.05             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 8.54              | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 8.38              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.21              | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 % | 6.70              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 6.28              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.75       | kN  | Min : | 1.69       | kN  | Punch NO.; | 27 |
| SD :                 | 1.90       | %   | Max : | 1.84       | kN  | Punch NO.; | 52 |
| Mean Main Pressure : | 7.99       | kN  | Min : | 6.81       | kN  | Punch NO.; | 36 |
| SD :                 | 6.20       | %   | Max : | 9.22       | kN  | Punch NO.; | 2  |
| Filling Depth :      | Increase ; | 348 | Times | Decrease ; | 354 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.55    | mm  | Pre Thickness ;  | 1.55              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.59    | mm  | Main Thickness ; | 1.59              | mm           |
| Filling Depth :               | 4.16    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.821   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 68      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 152     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.34 | kN |             |    |
|                  | Max ;     | 0.38 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                    |         |       |              |         |
|-----------------|--------------------|---------|-------|--------------|---------|
| Total Counter : | 6251218            | Tablets |       |              |         |
| A.W.C Counter : | 6217796            | Tablets |       |              |         |
| HEP :           | 2436 (0.0%)        | Tablets | LEP : | 64897 (1.0%) | Tablets |
| GOOD :          | 6150463<br>(98.9%) | Tablets |       |              |         |