

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP400MG STIU2B00
User ID :	os190259	Batch NO. :	IUM2500940
Print Interval :	50000 Tablets	Running Time :	09 Hour 24 Min 16 Sec
Total Counter :	1461452	Total Running Time :	2314 Hour 8 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	58	%	Manual ;	57	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	26.38	kN	No		
HEP :	20.0 %	25.33	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	21.53	kN	Yes 20	Times	
Ref :		21.11	kN			
LCP :	2.0 %	20.68	kN	Yes 20	Times	
LEP :	20.0 %	16.88	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	15.83	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.50	kN	Min :	2.32	kN	Punch NO.;	46
SD :	3.50	%	Max :	2.68	kN	Punch NO.;	41
Mean Main Pressure :	21.34	kN	Min :	19.39	kN	Punch NO.;	45
SD :	3.80	%	Max :	22.88	kN	Punch NO.;	41
Filling Depth :	Increase ; 49	Times	Decrease ; 67	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	139920	Tablets/hour
Feeder : Status ;	Auto		Speed ;	13	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.80	mm	Pre Thickness ;	2.81	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.40	mm
Filling Depth :	7.72	mm	Current Cam :	Cam C	
Main Air pressure :	0.840	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	15	Min			
Lower Punch Remain Time (S2):	27	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.07	kN		
	Max ;	0.18	kN	Max Punch ;	32
Lower Punch :	Average ;	0.05	kN		
	Max ;	0.15	kN	Max Punch ;	35
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.10	kN	Max Punch ;	35

Tablet Data (N/A)

Total Counter :	1461452	Tablets			
A.W.C Counter :	1437373	Tablets			
HEP :	1818 (0.1%)	Tablets	LEP :	707 (0.0%)	Tablets
GOOD :	1434848 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU BP400MG STIU2B00
User ID :	os190259	Batch NO. :	IUM2500940
Print Interval :	50000 Tablets	Running Time :	09 Hour 24 Min 17 Sec
Total Counter :	1461490	Total Running Time :	2314 Hour 8 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	62	%	Manual ;	50	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	27.07	kN	No		
HEP :	20.0 %	25.99	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	22.09	kN	Yes 20	Times	
Ref :		21.66	kN			
LCP :	2.0 %	21.23	kN	Yes 20	Times	
LEP :	20.0 %	17.33	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	16.24	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.72	kN	Min :	2.46	kN	Punch NO.;	2
SD :	4.90	%	Max :	3.08	kN	Punch NO.;	34
Mean Main Pressure :	21.69	kN	Min :	19.69	kN	Punch NO.;	8
SD :	6.00	%	Max :	25.37	kN	Punch NO.;	41
Filling Depth :	Increase ;	33	Times	Decrease ;	51	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	139920	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.80	mm	Pre Thickness ;	2.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.41	mm
Filling Depth :	8.21	mm	Current Cam :	Cam C	
Main Air pressure :	0.840	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	15	Min			
Lower Punch Remain Time (S2):	27	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.18	kN	Max Punch ;	41
Lower Punch :	Average ;	0.05	kN		
	Max ;	0.15	kN	Max Punch ;	35
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.10	kN	Max Punch ;	35

Tablet Data (N/A)

Total Counter :	1461490	Tablets			
A.W.C Counter :	1437373	Tablets			
HEP :	1818 (0.1%)	Tablets	LEP :	707 (0.0%)	Tablets
GOOD :	1434848 (99.8%)	Tablets			