

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 400 MG,STIQ2A00
User ID :	sj97904	Batch NO. :	IQM2502340
Print Interval :	50000 Tablets	Running Time :	07 Hour 24 Min 01 Sec
Total Counter :	700180	Total Running Time :	2357 Hour 1 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	143	%	Manual ;	20	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		810000	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	53.68	kN	No		
HEP :	20.0 %	51.53	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	43.80	kN	Yes 20	Times	
Ref :		42.94	kN			
LCP :	2.0 %	42.08	kN	Yes 20	Times	
LEP :	20.0 %	34.35	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	32.21	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	13.11	kN	Min :	11.48	kN	Punch NO.;	29
		5.00	%	Max :	14.75	kN	Punch NO.;	37
Mean Main Pressure :		43.03	kN	Min :	40.39	kN	Punch NO.;	29
	SD :	2.60	%	Max :	45.61	kN	Punch NO.;	37
Filling Depth :	Increase ;	9	Times	Decrease ;	33	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	21	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.00	mm	Pre Thickness ;	2.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.00	mm	Main Thickness ;	2.01	mm
Filling Depth :	6.85	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.828	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	3	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.16	kN		
	Max ;	0.27	kN	Max Punch ;	3
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.11	kN	Max Punch ;	24
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	700180	Tablets			
A.W.C Counter :	694645	Tablets			
HEP :	19 (0.0%)	Tablets	LEP :	119 (0.0%)	Tablets
GOOD :	694507 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 400 MG,STIQ2A00
User ID :	sj97904	Batch NO. :	IQM2502340
Print Interval :	50000 Tablets	Running Time :	07 Hour 24 Min 01 Sec
Total Counter :	700180	Total Running Time :	2357 Hour 1 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	150	%	Manual ;	21	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		810000	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	54.10	kN	No		
HEP :	20.0 %	51.93	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	44.14	kN	Yes 20	Times	
Ref :		43.28	kN			
LCP :	2.0 %	42.41	kN	Yes 20	Times	
LEP :	20.0 %	34.62	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	32.46	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	13.27	kN	Min :	11.32	kN	Punch NO.;	15
		10.10	%	Max :	17.00	kN	Punch NO.;	36
Mean Main Pressure :		43.57	kN	Min :	40.63	kN	Punch NO.;	14
	SD :	4.20	%	Max :	48.05	kN	Punch NO.;	36
Filling Depth :	Increase ;	46	Times	Decrease ;	10	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.00	mm	Pre Thickness ;	2.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.00	mm	Main Thickness ;	2.01	mm
Filling Depth :	7.12	mm	Current Cam :	Cam C	
Main Air pressure :	0.828	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	3	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.16	kN		
	Max ;	0.26	kN	Max Punch ;	1
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.11	kN	Max Punch ;	24
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	700180	Tablets			
A.W.C Counter :	694645	Tablets			
HEP :	19 (0.0%)	Tablets	LEP :	119 (0.0%)	Tablets
GOOD :	694507 (100.0%)	Tablets			