

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG,STIQ2A00   |
| User ID :        | sj97904       | Batch NO. :          | IQM2502340            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 31 Min 07 Sec |
| Total Counter :  | 805680        | Total Running Time : | 2358 Hour 7 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 152      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 810000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 53.68    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 51.53    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 43.80    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 42.94    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 42.08    | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 %      | 34.35    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 32.21    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 13.46 | kN    | Min :      | 11.46 | kN    | Punch NO.; | 14 |
|                      |            | 9.80  | %     | Max :      | 16.03 | kN    | Punch NO.; | 46 |
| Mean Main Pressure : |            | 43.47 | kN    | Min :      | 39.96 | kN    | Punch NO.; | 14 |
|                      | SD :       | 4.60  | %     | Max :      | 47.26 | kN    | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 14    | Times | Decrease ; | 3     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 23                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.00    | mm  | Main Thickness ; | 2.01              | mm           |
| Filling Depth :               | 7.33    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.824   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 20      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 20 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 805680             | Tablets |       |            |         |
| A.W.C Counter : | 798187             | Tablets |       |            |         |
| HEP :           | 32 (0.0%)          | Tablets | LEP : | 120 (0.0%) | Tablets |
| GOOD :          | 798035<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG,STIQ2A00   |
| User ID :        | sj97904       | Batch NO. :          | IQM2502340            |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 31 Min 07 Sec |
| Total Counter :  | 805680        | Total Running Time : | 2358 Hour 7 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 150      | %       | Manual ;           | 21           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 810000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 54.10    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 51.93    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 44.14    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 43.28    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 42.41    | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 %      | 34.62    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 32.46    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 12.65 | kN    | Min :      | 10.22 | kN    | Punch NO.; | 9  |
|                      |            | 11.50 | %     | Max :      | 17.25 | kN    | Punch NO.; | 46 |
| Mean Main Pressure : |            | 42.86 | kN    | Min :      | 39.62 | kN    | Punch NO.; | 9  |
|                      | SD :       | 4.40  | %     | Max :      | 48.30 | kN    | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 7     | Times | Decrease ; | 2     | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.00    | mm  | Main Thickness ; | 2.01              | mm           |
| Filling Depth :               | 7.10    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.824   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 20      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.53 | kN | Max Punch ; | 19 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 805680             | Tablets |       |            |         |
| A.W.C Counter : | 798187             | Tablets |       |            |         |
| HEP :           | 32 (0.0%)          | Tablets | LEP : | 120 (0.0%) | Tablets |
| GOOD :          | 798035<br>(100.0%) | Tablets |       |            |         |