

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Leveti1000mgSTLC2D00
User ID :	rp98927	Batch NO. :	LCM2502910
Print Interval :	50000 Tablets	Running Time :	04 Hour 34 Min 50 Sec
Total Counter :	450164	Total Running Time :	2376 Hour 34 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	170	%	Manual ;	20	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		660000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	48.12	kN	No		
HEP :	20.0 %	46.20	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	39.27	kN	Yes 20	Times	
Ref :		38.50	kN			
LCP :	2.0 %	37.73	kN	Yes 20	Times	
LEP :	20.0 %	30.80	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	28.87	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	2.13	kN	Min :	2.08	kN	Punch NO.;	30
		1.10	%	Max :	2.17	kN	Punch NO.;	1
Mean Main Pressure :	SD :	38.63	kN	Min :	36.88	kN	Punch NO.;	33
		2.00	%	Max :	40.21	kN	Punch NO.;	13
Filling Depth :	Increase ;	71	Times	Decrease ;	179	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	6.60	mm	Main Thickness ;	3.61	mm
Filling Depth :	10.13	mm	Current Cam :	Cam C	
Main Air pressure :	0.819	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	19	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.15	kN	Max Punch ;	6
Ejection Force :	Average ;	0.01	kN		
	Max ;	0.02	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	450164	Tablets			
A.W.C Counter :	256274	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	256274 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name : VANTIX P720D User ID : rp98927 Print Interval : 50000 Tablets Total Counter : 450164 Station : Station 2	Product Name : Leveti1000mgSTLC2D00 Batch NO. : LCM2502910 Running Time : 04 Hour 34 Min 50 Sec Total Running Time : 2376 Hour 34 Min
--	--

Setting Value

Feeder :	Auto ;	170	%	Manual ;	21	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
Target Quantity :		8.5-14mm			Max :	12.00 mm
Air Pressure Low Limit :		660000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	50.01	kN	No		
HEP :	20.0 %	48.01	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	40.81	kN	Yes 20	Times	
Ref :		40.01	kN			
LCP :	2.0 %	39.21	kN	Yes 20	Times	
LEP :	20.0 %	32.01	kN	Yes 15	Rotations 20	Tablets
LSP :	25.0 %	30.01	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.25	kN	Min :	2.19	kN	Punch NO.;	30
		1.70	%	Max :	2.36	kN	Punch NO.;	40
Mean Main Pressure :		40.08	kN	Min :	38.53	kN	Punch NO.;	17
	SD :	2.50	%	Max :	43.31	kN	Punch NO.;	13
Filling Depth :	Increase ;	61	Times	Decrease ;	97	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	6.75	mm	Main Thickness ;	3.76	mm
Filling Depth :	10.25	mm	Current Cam :	Cam C	
Main Air pressure :	0.819	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	19	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.15	kN	Max Punch ;	6
Ejection Force :	Average ;	0.01	kN		
	Max ;	0.02	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	450164	Tablets			
A.W.C Counter :	256274	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	256274 (100.0%)	Tablets			