

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Leveti1000mgSTLC2D00  |
| User ID :        | rp98927       | Batch NO. :          | LCM2503210            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 35 Min 39 Sec |
| Total Counter :  | 150034        | Total Running Time : | 2394 Hour 11 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 181      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 660000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 %      | 48.26    | kN      | No                 |              |           |
| HEP :                              | 20.0 %      | 46.33    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 39.38    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 38.60    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 37.83    | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 %      | 30.88    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 %      | 28.95    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.14           | kN    | Min :      | 2.09  | kN    | Punch NO.; | 23 |
| SD :                 | 1.20           | %     | Max :      | 2.21  | kN    | Punch NO.; | 42 |
| Mean Main Pressure : | 38.78          | kN    | Min :      | 37.09 | kN    | Punch NO.; | 33 |
| SD :                 | 1.80           | %     | Max :      | 40.11 | kN    | Punch NO.; | 11 |
| Filling Depth :      | Increase ; 106 | Times | Decrease ; | 122   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 27       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 9.00    | mm  | Pre Thickness ;  | 6.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.64    | mm  | Main Thickness ; | 3.65     | mm           |
| Filling Depth :               | 10.24   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.819   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 6  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 150034             | Tablets |       |        |         |
| A.W.C Counter : | 148066             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 148066<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Leveti1000mgSTLC2D00  |
| User ID :        | rp98927       | Batch NO. :          | LCM2503210            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 35 Min 39 Sec |
| Total Counter :  | 150034        | Total Running Time : | 2394 Hour 11 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 185       | %       | Manual ;           | 21           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 660000    | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80        | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30        | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 51.35     | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 49.30     | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 41.90     | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 41.08     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 40.26     | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 % | 32.87     | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 30.81     | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 40.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.33       | kN | Min : | 2.27       | kN  | Punch NO.; | 2  |
| SD :                 | 1.20       | %  | Max : | 2.41       | kN  | Punch NO.; | 40 |
| Mean Main Pressure : | 40.81      | kN | Min : | 39.30      | kN  | Punch NO.; | 15 |
| SD :                 | 2.10       | %  | Max : | 43.29      | kN  | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 84 | Times | Decrease ; | 127 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 9.00    | mm  | Pre Thickness ;  | 6.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.78    | mm  | Main Thickness ; | 3.79              | mm           |
| Filling Depth :               | 10.18   | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.819   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 5       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 6  |
| Ejection Force : | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 150034             | Tablets |       |        |         |
| A.W.C Counter : | 148066             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 148066<br>(100.0%) | Tablets |       |        |         |