

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA100MG STSH4B00
User ID :	sj97904	Batch NO. :	SHM2503810
Print Interval :	50000 Tablets	Running Time :	00 Hour 46 Min 13 Sec
Total Counter :	200166	Total Running Time :	2407 Hour 21 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	54	%	Manual ;	20	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	22.20	kN	No		
HEP :	20.0 %	21.32	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	18.12	kN	Yes 15	Times	
Ref :		17.76	kN			
LCP :	2.0 %	17.41	kN	Yes 15	Times	
LEP :	20.0 %	14.21	kN	Yes 10	Rotations 15	Tablets
LSP :	25.0 %	13.32	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.78	kN	Min :	1.73	kN	Punch NO.;	60
SD :	0.90	%	Max :	1.81	kN	Punch NO.;	10
Mean Main Pressure :	17.62	kN	Min :	16.39	kN	Punch NO.;	60
SD :	2.30	%	Max :	18.60	kN	Punch NO.;	75
Filling Depth :	Increase ; 6	Times	Decrease ; 4	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	258720	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.00	mm	Pre Thickness ;	4.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	5.54	mm	Main Thickness ;	2.54	mm
Filling Depth :	7.94	mm	Current Cam :	Cam C	
Main Air pressure :	0.837	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	4	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.46	kN		
	Max ;	0.64	kN	Max Punch ;	67
Lower Punch :	Average ;	0.17	kN		
	Max ;	0.30	kN	Max Punch ;	39
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200166	Tablets			
A.W.C Counter :	196224	Tablets			
HEP :	26 (0.0%)	Tablets	LEP :	482 (0.2%)	Tablets
GOOD :	195716 (99.7%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA100MG STSH4B00
User ID :	sj97904	Batch NO. :	SHM2503810
Print Interval :	50000 Tablets	Running Time :	00 Hour 46 Min 13 Sec
Total Counter :	200166	Total Running Time :	2407 Hour 21 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	50	%	Manual ;	21	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	25.0 %	23.01	kN	No		
HEP :	20.0 %	22.09	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	18.78	kN	Yes 15	Times	
Ref :		18.41	kN			
LCP :	2.0 %	18.04	kN	Yes 15	Times	
LEP :	20.0 %	14.73	kN	Yes 10	Rotations 15	Tablets
LSP :	25.0 %	13.81	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.97	kN	Min :	1.92	kN	Punch NO.;	24
SD :	1.30	%	Max :	2.05	kN	Punch NO.;	40
Mean Main Pressure :	18.40	kN	Min :	16.90	kN	Punch NO.;	18
SD :	3.30	%	Max :	20.02	kN	Punch NO.;	57
Filling Depth :	Increase ;	8	Times	Decrease ;	10	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	258720	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.00	mm	Pre Thickness ;	4.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	5.70	mm	Main Thickness ;	2.70	mm
Filling Depth :	7.88	mm	Current Cam :	Cam C	
Main Air pressure :	0.837	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	4	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.46	kN		
	Max ;	0.64	kN	Max Punch ;	67
Lower Punch :	Average ;	0.17	kN		
	Max ;	0.30	kN	Max Punch ;	39
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200166	Tablets			
A.W.C Counter :	196224	Tablets			
HEP :	26 (0.0%)	Tablets	LEP :	482 (0.2%)	Tablets
GOOD :	195716 (99.7%)	Tablets			