

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Allo200mg STUB3D00    |
| User ID :        | cb83530       | Batch NO. :          | UBM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 31 Min 40 Sec |
| Total Counter :  | 550174        | Total Running Time : | 2518 Hour 45 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 150      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 80       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 28.0 %      | 13.94    | kN      | No                 |              |           |
| HEP :                              | 23.0 %      | 13.40    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 11.11    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.89    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.68    | kN      | Yes 15             | Times        |           |
| LEP :                              | 23.0 %      | 8.39     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 28.0 %      | 7.84     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.54           | kN    | Min :      | 1.52  | kN    | Punch NO.; | 36 |
| SD :                 | 0.60           | %     | Max :      | 1.56  | kN    | Punch NO.; | 6  |
| Mean Main Pressure : | 11.47          | kN    | Min :      | 10.50 | kN    | Punch NO.; | 19 |
| SD :                 | 5.80           | %     | Max :      | 13.49 | kN    | Punch NO.; | 37 |
| Filling Depth :      | Increase ; 727 | Times | Decrease ; | 711   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 30                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.00              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.08    | mm  | Main Thickness ; | 3.08              | mm           |
| Filling Depth :               | 7.73    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.733   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 32      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 62      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 5  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 550174            | Tablets |       |            |         |
| A.W.C Counter : | 549027            | Tablets |       |            |         |
| HEP :           | 718 (0.1%)        | Tablets | LEP : | 106 (0.0%) | Tablets |
| GOOD :          | 548203<br>(99.8%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Allo200mg STUB3D00    |
| User ID :        | cb83530       | Batch NO. :          | UBM2500320            |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 31 Min 40 Sec |
| Total Counter :  | 550174        | Total Running Time : | 2518 Hour 45 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 138      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 80       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 28.0 %      | 12.83    | kN      | No                 |              |           |
| HEP :                              | 23.0 %      | 12.33    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 10.23    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.03    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 9.83     | kN      | Yes 15             | Times        |           |
| LEP :                              | 23.0 %      | 7.72     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 28.0 %      | 7.22     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.18           | kN    | Min :      | 1.16  | kN    | Punch NO.; | 33 |
| SD :                 | 0.90           | %     | Max :      | 1.20  | kN    | Punch NO.; | 34 |
| Mean Main Pressure : | 9.85           | kN    | Min :      | 8.60  | kN    | Punch NO.; | 52 |
| SD :                 | 5.60           | %     | Max :      | 11.12 | kN    | Punch NO.; | 17 |
| Filling Depth :      | Increase ; 699 | Times | Decrease ; | 697   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.35    | mm  | Main Thickness ; | 3.35     | mm           |
| Filling Depth :               | 8.30    | mm  | Current Cam :    | Cam C    |              |
|                               |         |     |                  | 8.5-14mm |              |
| Main Air pressure :           | 0.733   | MPa |                  |          |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 32      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 62      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 5  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 550174            | Tablets |       |            |         |
| A.W.C Counter : | 549027            | Tablets |       |            |         |
| HEP :           | 718 (0.1%)        | Tablets | LEP : | 106 (0.0%) | Tablets |
| GOOD :          | 548203<br>(99.8%) | Tablets |       |            |         |