

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA<br>50MG,STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2503150              |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 50 Min 16 Sec   |
| Total Counter :  | 5907390       | Total Running Time : | 2549 Hour 31 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 55                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 50                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 28.0 %      | 19.00             | kN      | No                 |              |           |
| HEP :                              | 23.0 %      | 18.26             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 15.14             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.85             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.55             | kN      | Yes 15             | Times        |           |
| LEP :                              | 23.0 %      | 11.43             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 28.0 %      | 10.69             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.79          | kN    | Min :      | 1.75  | kN    | Punch NO.; | 2  |
| SD :                 | 0.90          | %     | Max :      | 1.83  | kN    | Punch NO.; | 65 |
| Mean Main Pressure : | 14.97         | kN    | Min :      | 14.08 | kN    | Punch NO.; | 41 |
| SD :                 | 3.50          | %     | Max :      | 16.36 | kN    | Punch NO.; | 69 |
| Filling Depth :      | Increase ; 51 | Times | Decrease ; | 154   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.03    | mm  | Main Thickness ; | 2.03              | mm           |
| Filling Depth :               | 6.77    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.823   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 29      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.39 | kN |             |    |
|                  | Max ;     | 0.63 | kN | Max Punch ; | 34 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 58 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 5907390             | Tablets |       |          |         |
| A.W.C Counter : | 5898101             | Tablets |       |          |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 4 (0.0%) | Tablets |
| GOOD :          | 5898097<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA<br>50MG,STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2503150              |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 50 Min 16 Sec   |
| Total Counter :  | 5907390       | Total Running Time : | 2549 Hour 31 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 55                | %       | Manual ;           | 40           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 50                | Min     | Supply Time ;      | 1.0          | Sec        |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 28.0 %      | 19.04             | kN      | No                 |              |            |
| HEP :                              | 23.0 %      | 18.30             | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.0 %       | 15.18             | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 14.88             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 14.58             | kN      | Yes                | 15 Times     |            |
| LEP :                              | 23.0 %      | 11.46             | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 28.0 %      | 10.71             | kN      | No                 |              |            |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.40  | kN    | Min :      | 1.37  | kN    | Punch NO.; | 50 |
|                      |            | 1.40  | %     | Max :      | 1.46  | kN    | Punch NO.; | 72 |
| Mean Main Pressure : | SD :       | 12.96 | kN    | Min :      | 11.18 | kN    | Punch NO.; | 36 |
|                      |            | 6.00  | %     | Max :      | 14.88 | kN    | Punch NO.; | 73 |
| Filling Depth :      | Increase ; | 48    | Times | Decrease ; | 59    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.06    | mm  | Main Thickness ; | 2.06              | mm           |
| Filling Depth :               | 6.90    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.823   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 29      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.38 | kN |             |    |
|                  | Max ;     | 0.70 | kN | Max Punch ; | 34 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 58 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 5907390             | Tablets |       |          |         |
| A.W.C Counter : | 5898101             | Tablets |       |          |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 4 (0.0%) | Tablets |
| GOOD :          | 5898097<br>(100.0%) | Tablets |       |          |         |