

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU BP600MG<br>STIU2C00 |
| User ID :        | rk99163       | Batch NO. :          | IUM2504460              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 33 Min 07 Sec   |
| Total Counter :  | 1000238       | Total Running Time : | 2678 Hour 33 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 250               | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 1350000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 5.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.99             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 28.84             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 23.53             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 23.07             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 22.61             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 17.30             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 16.15             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.90       | kN  | Min : | 3.87       | kN  | Punch NO.; | 41 |
| SD :                 | 0.50       | %   | Max : | 3.95       | kN  | Punch NO.; | 4  |
| Mean Main Pressure : | 23.00      | kN  | Min : | 21.54      | kN  | Punch NO.; | 51 |
| SD :                 | 3.40       | %   | Max : | 24.63      | kN  | Punch NO.; | 16 |
| Filling Depth :      | Increase ; | 236 | Times | Decrease ; | 184 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 38       | RPM          |
| Pre Upper Punch Position ;    | 3.98    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.60    | mm  | Pre Thickness ;  | 4.62     | mm           |
| Main Upper Punch Position ;   | 3.98    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 8.42    | mm  | Main Thickness ; | 4.44     | mm           |
| Filling Depth :               | 11.55   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 4       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 19 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 5  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1000238            | Tablets |       |            |         |
| A.W.C Counter : | 986628             | Tablets |       |            |         |
| HEP :           | 86 (0.0%)          | Tablets | LEP : | 315 (0.0%) | Tablets |
| GOOD :          | 986227<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU BP600MG<br>STIU2C00 |
| User ID :        | rk99163       | Batch NO. :          | IUM2504460              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 33 Min 07 Sec   |
| Total Counter :  | 1000238       | Total Running Time : | 2678 Hour 33 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 220               | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 1350000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 5.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 29.97             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 28.81             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 23.51             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 23.05             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 22.59             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 17.29             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 16.14             | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.94       | kN  | Min : | 3.90       | kN  | Punch NO.; | 34 |
| SD :                 | 0.60       | %   | Max : | 4.03       | kN  | Punch NO.; | 48 |
| Mean Main Pressure : | 24.08      | kN  | Min : | 20.63      | kN  | Punch NO.; | 45 |
| SD :                 | 7.00       | %   | Max : | 27.42      | kN  | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 212 | Times | Decrease ; | 151 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 33       | RPM          |
| Pre Upper Punch Position ;    | 4.01    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.60    | mm  | Pre Thickness ;  | 4.59     | mm           |
| Main Upper Punch Position ;   | 3.98    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 8.42    | mm  | Main Thickness ; | 4.44     | mm           |
| Filling Depth :               | 10.94   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.816   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 4       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.10 | kN |             |   |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 4 |
| Lower Punch :    | Average ; | 0.07 | kN |             |   |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 5 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1000238            | Tablets |       |            |         |
| A.W.C Counter : | 986628             | Tablets |       |            |         |
| HEP :           | 86 (0.0%)          | Tablets | LEP : | 315 (0.0%) | Tablets |
| GOOD :          | 986227<br>(100.0%) | Tablets |       |            |         |