

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Gaba 600mg STKS1A00   |
| User ID :        | rp98927       | Batch NO. :          | KSM2500160            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 28 Min 11 Sec |
| Total Counter :  | 50056         | Total Running Time : | 2829 Hour 29 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |                        |         |                    |              |           |
|------------------------------------|-------------|------------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 120                    | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min :        | 10.50 mm  |
| Target Quantity :                  |             | 1460000                | Tablets |                    | Max :        | 16.00 mm  |
| Air Pressure Low Limit :           |             | 0.000                  | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                        |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50                     | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60                     | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0                   | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                      | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00                   | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00                   | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90                   | kN      |                    |              |           |
|                                    |             |                        |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.01                  | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.36                  | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 13.35                  | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.08                  | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.82                  | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.81                   | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 9.16                   | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                        |         | Yes                |              |           |
| Pre HSP                            |             | 40.00                  | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                        |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                        |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 3.50          | kN    | Min :         | 3.40  | kN | Punch NO.; | 36 |
| SD :                 | 1.50          | %     | Max :         | 3.65  | kN | Punch NO.; | 13 |
| Mean Main Pressure : | 13.18         | kN    | Min :         | 12.48 | kN | Punch NO.; | 29 |
| SD :                 | 2.40          | %     | Max :         | 13.99 | kN | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 11 | Times | Decrease ; 19 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20                     | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.20    | mm  | Pre Thickness ;  | 5.20                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.30    | mm  | Main Thickness ; | 4.30                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 11.15   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.821   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 28      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 47      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 50056             | Tablets |       |        |         |
| A.W.C Counter : | 46658             | Tablets |       |        |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 46658<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Gaba 600mg STKS1A00   |
| User ID :        | rp98927       | Batch NO. :          | KSM2500160            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 28 Min 11 Sec |
| Total Counter :  | 50056         | Total Running Time : | 2829 Hour 29 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |                        |         |                    |              |           |
|------------------------------------|--------|------------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 120                    | %       | Manual ;           | 27           | RPM       |
| Filling Cam :                      |        | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min :        | 10.50 mm  |
| Target Quantity :                  |        | 1460000                | Tablets |                    | Max :        | 16.00 mm  |
| Air Pressure Low Limit :           |        | 0.000                  | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                        |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                     | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60                     | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                      | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00                   | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00                   | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90                   | kN      |                    |              |           |
|                                    |        |                        |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 18.29                  | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 17.59                  | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 14.35                  | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 14.07                  | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 13.79                  | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 10.55                  | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 9.85                   | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                        |         | Yes                |              |           |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                        |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 3.01         | kN    | Min :        | 2.86  | kN | Punch NO.; | 36 |
| SD :                 | 2.10         | %     | Max :        | 3.17  | kN | Punch NO.; | 19 |
| Mean Main Pressure : | 13.96        | kN    | Min :        | 13.20 | kN | Punch NO.; | 4  |
| SD :                 | 3.00         | %     | Max :        | 15.15 | kN | Punch NO.; | 27 |
| Filling Depth :      | Increase ; 6 | Times | Decrease ; 6 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20                     | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.20    | mm  | Pre Thickness ;  | 5.20                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.48    | mm  | Main Thickness ; | 4.48                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 11.26   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.821   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 28      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 47      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 50056             | Tablets |       |        |         |
| A.W.C Counter : | 46658             | Tablets |       |        |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 46658<br>(100.0%) | Tablets |       |        |         |