

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | sj97904       | Batch NO. :          | KQM2509940            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 27 Min 47 Sec |
| Total Counter :  | 1375216       | Total Running Time : | 2900 Hour 50 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 52       | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 38.39    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 36.91    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 30.12    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 29.53    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 28.94    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.15    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 20.67    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.92       | kN  | Min : | 3.59       | kN  | Punch NO.; | 51 |
| SD :                 | 4.00       | %   | Max : | 4.34       | kN  | Punch NO.; | 4  |
| Mean Main Pressure : | 29.17      | kN  | Min : | 23.60      | kN  | Punch NO.; | 44 |
| SD :                 | 5.60       | %   | Max : | 33.66      | kN  | Punch NO.; | 34 |
| Filling Depth :      | Increase ; | 423 | Times | Decrease ; | 424 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.80    | mm  | Main Thickness ; | 0.80     | mm           |
| Filling Depth :               | 3.19    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.805   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.20 | kN |             |   |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 1 |
| Lower Punch :    | Average ; | 0.14 | kN |             |   |
|                  | Max ;     | 0.40 | kN | Max Punch ; | 4 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1375216             | Tablets |       |           |         |
| A.W.C Counter : | 1368889             | Tablets |       |           |         |
| HEP :           | 22 (0.0%)           | Tablets | LEP : | 63 (0.0%) | Tablets |
| GOOD :          | 1368804<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | sj97904       | Batch NO. :          | KQM2509940            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 27 Min 48 Sec |
| Total Counter :  | 1375266       | Total Running Time : | 2900 Hour 50 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 52       | %       | Manual ;           | 27           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 37.77    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 36.32    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 29.64    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 29.06    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 28.48    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 21.79    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 20.34    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.93  | kN    | Min :      | 3.61  | kN    | Punch NO.; | 28 |
|                      |            | 3.30  | %     | Max :      | 4.24  | kN    | Punch NO.; | 14 |
| Mean Main Pressure : |            | 29.15 | kN    | Min :      | 25.18 | kN    | Punch NO.; | 44 |
|                      | SD :       | 3.90  | %     | Max :      | 31.58 | kN    | Punch NO.; | 34 |
| Filling Depth :      | Increase ; | 248   | Times | Decrease ; | 254   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.84    | mm  | Main Thickness ; | 0.85              | mm           |
| Filling Depth :               | 3.70    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.806   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 38      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 17 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.40 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1375266             | Tablets |       |           |         |
| A.W.C Counter : | 1368889             | Tablets |       |           |         |
| HEP :           | 22 (0.0%)           | Tablets | LEP : | 63 (0.0%) | Tablets |
| GOOD :          | 1368804<br>(100.0%) | Tablets |       |           |         |