

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Panto 40 mg STKQ1B00  |
| User ID :        | am98755       | Batch NO. :          | KQM2510040            |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 46 Min 06 Sec |
| Total Counter :  | 2182190       | Total Running Time : | 2918 Hour 43 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 42       | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 10       | Min     | Supply Time ;      | 4.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 38.06    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 36.60    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 29.86    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 29.28    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 28.69    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 21.96    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 20.50    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 4.12          | kN    | Min :         | 3.50  | kN | Punch NO.; | 5  |
| SD :                 | 4.80          | %     | Max :         | 4.56  | kN | Punch NO.; | 26 |
| Mean Main Pressure : | 29.43         | kN    | Min :         | 22.65 | kN | Punch NO.; | 44 |
| SD :                 | 6.80          | %     | Max :         | 34.35 | kN | Punch NO.; | 76 |
| Filling Depth :      | Increase ; 97 | Times | Decrease ; 77 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 203280   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.81    | mm  | Main Thickness ; | 0.82     | mm           |
| Filling Depth :               | 3.88    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.714   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.24 | kN |             |    |
|                  | Max ;     | 0.33 | kN | Max Punch ; | 50 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.46 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 2182190            | Tablets |       |             |         |
| A.W.C Counter : | 2168529            | Tablets |       |             |         |
| HEP :           | 456 (0.0%)         | Tablets | LEP : | 2382 (0.1%) | Tablets |
| GOOD :          | 2165691<br>(99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                                                                                                                                                                         |                                                                                                                                                                            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Machine Name :</b> VANTIX P720D<br><b>User ID :</b> am98755<br><b>Print Interval :</b> 50000 Tablets<br><b>Total Counter :</b> 2182190<br><b>Station :</b> Station 2 | <b>Product Name :</b> Panto 40 mg STKQ1B00<br><b>Batch NO. :</b> KQM2510040<br><b>Running Time :</b> 10 Hour 46 Min 06 Sec<br><b>Total Running Time :</b> 2918 Hour 43 Min |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 42       | %       | Manual ;           | 25           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 10       | Min     | Supply Time ;      | 4.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 35.61    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 34.24    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 27.94    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 27.39    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 26.84    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 20.54    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 19.17    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.99  | kN    | Min :      | 3.39  | kN    | Punch NO.; | 3  |
|                      |            | 5.60  | %     | Max :      | 4.44  | kN    | Punch NO.; | 57 |
| Mean Main Pressure : |            | 27.29 | kN    | Min :      | 18.62 | kN    | Punch NO.; | 44 |
|                      | SD :       | 6.90  | %     | Max :      | 31.03 | kN    | Punch NO.; | 19 |
| Filling Depth :      | Increase ; | 32    | Times | Decrease ; | 28    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 203280   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.86    | mm  | Main Thickness ; | 0.87     | mm           |
| Filling Depth :               | 3.78    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.714   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.46 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 2182190            | Tablets |       |             |         |
| A.W.C Counter : | 2168529            | Tablets |       |             |         |
| HEP :           | 456 (0.0%)         | Tablets | LEP : | 2382 (0.1%) | Tablets |
| GOOD :          | 2165691<br>(99.9%) | Tablets |       |             |         |