

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Allo 200mg STUB3D00
User ID :	am98755	Batch NO. :	UBM2500520
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 32 Sec
Total Counter :	50004	Total Running Time :	2934 Hour 9 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	75	%	Manual ;	23	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
Target Quantity :		8.5-14mm			Max :	12.00 mm
Air Pressure Low Limit :		1500000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	20.36	kN	No		
HEP :	25.0 %	19.58	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	16.05	kN	Yes 15	Times	
Ref :		15.66	kN			
LCP :	2.5 %	15.27	kN	Yes 15	Times	
LEP :	25.0 %	11.75	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	10.96	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.21	kN	Min :	2.18	kN	Punch NO.;	48
SD :	0.30	%	Max :	2.22	kN	Punch NO.;	34
Mean Main Pressure :	15.81	kN	Min :	14.53	kN	Punch NO.;	33
SD :	3.30	%	Max :	16.76	kN	Punch NO.;	9
Filling Depth :	Increase ;	256	Times	Decrease ;	266	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.30	mm	Main Thickness ;	3.30	mm
Filling Depth :	7.92	mm	Current Cam :	Cam C	
Main Air pressure :	0.747	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	20	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.11	kN	Max Punch ;	32
Lower Punch :	Average ;	0.12	kN		
	Max ;	0.16	kN	Max Punch ;	44
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50004	Tablets			
A.W.C Counter :	46920	Tablets			
HEP :	156 (0.3%)	Tablets	LEP :	76 (0.2%)	Tablets
GOOD :	46688 (99.5%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Allo 200mg STUB3D00
User ID :	am98755	Batch NO. :	UBM2500520
Print Interval :	50000 Tablets	Running Time :	00 Hour 31 Min 32 Sec
Total Counter :	50004	Total Running Time :	2934 Hour 9 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	75	%	Manual ;	25	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		1500000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	19.21	kN	No		
HEP :	25.0 %	18.47	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	15.15	kN	Yes 15	Times	
Ref :		14.78	kN			
LCP :	2.5 %	14.41	kN	Yes 15	Times	
LEP :	25.0 %	11.08	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	10.34	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.15	kN	Min :	2.13	kN	Punch NO.;	25
		0.50	%	Max :	2.18	kN	Punch NO.;	53
Mean Main Pressure :		14.87	kN	Min :	13.05	kN	Punch NO.;	49
	SD :	5.30	%	Max :	16.52	kN	Punch NO.;	15
Filling Depth :	Increase ;	591	Times	Decrease ;	520	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	11	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.53	mm	Main Thickness ;	3.53	mm
Filling Depth :	8.61	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.747	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	20	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.08	kN		
	Max ;	0.11	kN	Max Punch ;	32
Lower Punch :	Average ;	0.12	kN		
	Max ;	0.16	kN	Max Punch ;	44
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50004	Tablets			
A.W.C Counter :	46920	Tablets			
HEP :	156 (0.3%)	Tablets	LEP :	76 (0.2%)	Tablets
GOOD :	46688 (99.5%)	Tablets			