

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Allo 200mg STUB3D00
User ID :	pb97604	Batch NO. :	UBM2500520
Print Interval :	50000 Tablets	Running Time :	01 Hour 42 Min 35 Sec
Total Counter :	200084	Total Running Time :	2935 Hour 20 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	140	%	Manual ;	23	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
Target Quantity :		8.5-14mm			Max :	12.00 mm
Air Pressure Low Limit :		1500000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	20.93	kN	No		
HEP :	25.0 %	20.13	kN	Yes	10	Rotations 15 Tablets
HCP :	2.5 %	16.50	kN	Yes	15	Times
Ref :		16.10	kN			
LCP :	2.5 %	15.70	kN	Yes	15	Times
LEP :	25.0 %	12.08	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	11.27	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.18	kN	Min :	2.16	kN	Punch NO.;	12
SD :	0.30	%	Max :	2.21	kN	Punch NO.;	1
Mean Main Pressure :	15.88	kN	Min :	12.44	kN	Punch NO.;	34
SD :	8.60	%	Max :	19.81	kN	Punch NO.;	21
Filling Depth :	Increase ;	331	Times	Decrease ;	340	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	28	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.27	mm	Main Thickness ;	3.27	mm
Filling Depth :	8.37	mm	Current Cam :	Cam C	
Main Air pressure :	0.754	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.22	kN		
	Max ;	0.26	kN	Max Punch ;	17
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.22	kN	Max Punch ;	38
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200084	Tablets			
A.W.C Counter :	191680	Tablets			
HEP :	699 (0.4%)	Tablets	LEP :	173 (0.1%)	Tablets
GOOD :	190808 (99.5%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name : VANTIX P720D User ID : pb97604 Print Interval : 50000 Tablets Total Counter : 200084 Station : Station 2	Product Name : Allo 200mg STUB3D00 Batch NO. : UBM2500520 Running Time : 01 Hour 42 Min 35 Sec Total Running Time : 2935 Hour 20 Min
--	---

Setting Value

Feeder :	Auto ;	128	%	Manual ;	25	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
Target Quantity :		8.5-14mm			Max :	12.00 mm
Air Pressure Low Limit :		1500000	Tablets			
Oil Lubrication :		0.400	MPa			
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	19.24	kN	No		
HEP :	25.0 %	18.50	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	15.17	kN	Yes 15	Times	
Ref :		14.80	kN			
LCP :	2.5 %	14.43	kN	Yes 15	Times	
LEP :	25.0 %	11.10	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	10.36	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.16	kN	Min :	2.12	kN	Punch NO.;	40
		0.80	%	Max :	2.21	kN	Punch NO.;	20
Mean Main Pressure :		14.66	kN	Min :	13.44	kN	Punch NO.;	51
	SD :	4.40	%	Max :	16.98	kN	Punch NO.;	32
Filling Depth :	Increase ;	694	Times	Decrease ;	592	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	7.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	6.53	mm	Main Thickness ;	3.53	mm
Filling Depth :	9.00	mm	Current Cam :	Cam C	
				8.5-14mm	
Main Air pressure :	0.754	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.22	kN		
	Max ;	0.26	kN	Max Punch ;	17
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.22	kN	Max Punch ;	38
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	200084	Tablets			
A.W.C Counter :	191680	Tablets			
HEP :	699 (0.4%)	Tablets	LEP :	173 (0.1%)	Tablets
GOOD :	190808 (99.5%)	Tablets			