

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	rp98927	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	01 Hour 30 Min 23 Sec
Total Counter :	300212	Total Running Time :	2947 Hour 42 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	40	%	Manual ;	23	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.29	kN	No		
HEP :	25.0 %	38.74	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	31.77	kN	Yes 15	Times	
Ref :		30.99	kN			
LCP :	2.5 %	30.22	kN	Yes 15	Times	
LEP :	25.0 %	23.25	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.70	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	4.07	kN	Min :	3.41	kN	Punch NO.;	17
		7.30	%	Max :	5.86	kN	Punch NO.;	70
Mean Main Pressure :	SD :	31.21	kN	Min :	26.89	kN	Punch NO.;	9
		6.50	%	Max :	43.66	kN	Punch NO.;	70
Filling Depth :	Increase ;	36	Times	Decrease ;	34	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	9	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.84	mm	Main Thickness ;	0.84	mm
Filling Depth :	3.82	mm	Current Cam :	Cam B	
Main Air pressure :	0.760	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	16	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.18	kN	Max Punch ;	18
Lower Punch :	Average ;	0.42	kN		
	Max ;	0.99	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	300212	Tablets			
A.W.C Counter :	287668	Tablets			
HEP :	36 (0.0%)	Tablets	LEP :	547 (0.2%)	Tablets
GOOD :	287085 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	rp98927	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	01 Hour 30 Min 23 Sec
Total Counter :	300212	Total Running Time :	2947 Hour 42 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	40	%	Manual ;	25	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	39.51	kN	No		
HEP :	25.0 %	37.99	kN	Yes	10	Rotations 15 Tablets
HCP :	2.5 %	31.15	kN	Yes	15	Times
Ref :		30.39	kN			
LCP :	2.5 %	29.63	kN	Yes	15	Times
LEP :	25.0 %	22.79	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	21.27	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.82	kN	Min :	2.26	kN	Punch NO.;	70
SD :	6.70	%	Max :	4.53	kN	Punch NO.;	46
Mean Main Pressure :	29.65	kN	Min :	21.54	kN	Punch NO.;	70
SD :	4.80	%	Max :	33.09	kN	Punch NO.;	46
Filling Depth :	Increase ;	6	Times	Decrease ;	22	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	9	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.77	mm	Main Thickness ;	0.77	mm
Filling Depth :	3.64	mm	Current Cam :	Cam B	
Main Air pressure :	0.760	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	16	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.18	kN	Max Punch ;	18
Lower Punch :	Average ;	0.42	kN		
	Max ;	0.96	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	300212	Tablets			
A.W.C Counter :	287668	Tablets			
HEP :	36 (0.0%)	Tablets	LEP :	547 (0.2%)	Tablets
GOOD :	287085 (99.8%)	Tablets			