

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	rp98927	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	01 Hour 45 Min 22 Sec
Total Counter :	350214	Total Running Time :	2947 Hour 57 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	46	%	Manual ;	23	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.29	kN	No		
HEP :	25.0 %	38.74	kN	Yes 10	Rotations 15	Tablets
HCP :	2.5 %	31.77	kN	Yes 15	Times	
Ref :		30.99	kN			
LCP :	2.5 %	30.22	kN	Yes 15	Times	
LEP :	25.0 %	23.25	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.70	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.00	kN	Min :	3.56	kN	Punch NO.;	55
SD :	3.90	%	Max :	4.39	kN	Punch NO.;	6
Mean Main Pressure :	30.79	kN	Min :	25.48	kN	Punch NO.;	52
SD :	4.00	%	Max :	33.24	kN	Punch NO.;	43
Filling Depth :	Increase ; 38	Times	Decrease ; 36	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	10	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.84	mm	Main Thickness ;	0.84	mm
Filling Depth :	4.15	mm	Current Cam :	Cam B	
Main Air pressure :	0.755	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.24	kN	Max Punch ;	30
Lower Punch :	Average ;	0.43	kN		
	Max ;	0.91	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	350214	Tablets			
A.W.C Counter :	336741	Tablets			
HEP :	46 (0.0%)	Tablets	LEP :	551 (0.2%)	Tablets
GOOD :	336144 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	rp98927	Batch NO. :	KQM2510240
Print Interval :	50000 Tablets	Running Time :	01 Hour 45 Min 22 Sec
Total Counter :	350214	Total Running Time :	2947 Hour 57 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	40	%	Manual ;	25	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		20	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	39.93	kN	No		
HEP :	25.0 %	38.40	kN	Yes	10	Rotations 15 Tablets
HCP :	2.5 %	31.49	kN	Yes	15	Times
Ref :		30.72	kN			
LCP :	2.5 %	29.95	kN	Yes	15	Times
LEP :	25.0 %	23.04	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	21.50	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.92	kN	Min :	3.57	kN	Punch NO.;	52
SD :	4.90	%	Max :	4.63	kN	Punch NO.;	46
Mean Main Pressure :	30.28	kN	Min :	25.61	kN	Punch NO.;	52
SD :	4.10	%	Max :	33.66	kN	Punch NO.;	46
Filling Depth :	Increase ;	8	Times	Decrease ;	27	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	9	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.77	mm	Main Thickness ;	0.77	mm
Filling Depth :	3.69	mm	Current Cam :	Cam B	
Main Air pressure :	0.755	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	16	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.24	kN	Max Punch ;	30
Lower Punch :	Average ;	0.43	kN		
	Max ;	0.91	kN	Max Punch ;	9
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	350214	Tablets			
A.W.C Counter :	336741	Tablets			
HEP :	46 (0.0%)	Tablets	LEP :	551 (0.2%)	Tablets
GOOD :	336144 (99.8%)	Tablets			