

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2510240              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 33 Min 23 Sec   |
| Total Counter :  | 500254        | Total Running Time : | 2948 Hour 45 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 46                | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 41.10             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 39.52             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 32.41             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 31.62             | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 30.83             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 23.71             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 22.13             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.05          | kN    | Min :      | 3.71  | kN    | Punch NO.; | 2  |
| SD :                 | 5.10          | %     | Max :      | 4.76  | kN    | Punch NO.; | 46 |
| Mean Main Pressure : | 32.01         | kN    | Min :      | 27.68 | kN    | Punch NO.; | 52 |
| SD :                 | 4.30          | %     | Max :      | 35.72 | kN    | Punch NO.; | 46 |
| Filling Depth :      | Increase ; 39 | Times | Decrease ; | 54    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.84    | mm  | Main Thickness ; | 0.84     | mm           |
| Filling Depth :               | 3.70    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.753   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.10 | kN |             |   |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 5 |
| Lower Punch :    | Average ; | 0.41 | kN |             |   |
|                  | Max ;     | 0.93 | kN | Max Punch ; | 9 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500254            | Tablets |       |            |         |
| A.W.C Counter : | 484332            | Tablets |       |            |         |
| HEP :           | 71 (0.0%)         | Tablets | LEP : | 573 (0.1%) | Tablets |
| GOOD :          | 483688<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2510240              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 33 Min 23 Sec   |
| Total Counter :  | 500254        | Total Running Time : | 2948 Hour 45 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 45                | %       | Manual ;           | 25           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec        |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 41.70             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 40.10             | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.5 %       | 32.88             | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 32.08             | kN      |                    |              |            |
| LCP :                              | 2.5 %       | 31.28             | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 24.06             | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 30.0 %      | 22.46             | kN      | No                 |              |            |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 1 Rotations |                   |         | Calculation Pass : | 1            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 4.07       | kN | Min : | 3.55       | kN | Punch NO.; | 58 |
| SD :                 | 4.80       | %  | Max : | 4.52       | kN | Punch NO.; | 53 |
| Mean Main Pressure : | 32.05      | kN | Min : | 27.21      | kN | Punch NO.; | 52 |
| SD :                 | 3.70       | %  | Max : | 34.57      | kN | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 11 | Times | Decrease ; | 30 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.77    | mm  | Main Thickness ; | 0.77     | mm           |
| Filling Depth :               | 3.79    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.753   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |   |
|------------------|-----------|------|----|-------------|---|
| Upper Punch :    | Average ; | 0.10 | kN |             |   |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 5 |
| Lower Punch :    | Average ; | 0.41 | kN |             |   |
|                  | Max ;     | 0.96 | kN | Max Punch ; | 9 |
| Ejection Force : | Average ; | 0.00 | kN |             |   |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500254            | Tablets |       |            |         |
| A.W.C Counter : | 484332            | Tablets |       |            |         |
| HEP :           | 71 (0.0%)         | Tablets | LEP : | 573 (0.1%) | Tablets |
| GOOD :          | 483688<br>(99.9%) | Tablets |       |            |         |