

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2510240              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 22 Min 01 Sec   |
| Total Counter :  | 1750882       | Total Running Time : | 2955 Hour 34 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 50                | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 39.76             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 38.23             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 31.19             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 30.58             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 29.97             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.94             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 21.41             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.83       | kN  | Min : | 3.38       | kN  | Punch NO.; | 22 |
| SD :                 | 4.20       | %   | Max : | 4.29       | kN  | Punch NO.; | 61 |
| Mean Main Pressure : | 29.70      | kN  | Min : | 26.27      | kN  | Punch NO.; | 47 |
| SD :                 | 4.00       | %   | Max : | 33.08      | kN  | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 493 | Times | Decrease ; | 479 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10                | RPM          |
| Pre Upper Punch Position ;    | 2.69    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.31              | mm           |
| Main Upper Punch Position ;   | 2.69    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.91              | mm           |
| Filling Depth :               | 4.03    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.802   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 65 |
| Lower Punch :    | Average ; | 0.48 | kN |             |    |
|                  | Max ;     | 0.92 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1750882            | Tablets |       |             |         |
| A.W.C Counter : | 1719061            | Tablets |       |             |         |
| HEP :           | 249 (0.0%)         | Tablets | LEP : | 1656 (0.1%) | Tablets |
| GOOD :          | 1717156<br>(99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2510240              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 22 Min 01 Sec   |
| Total Counter :  | 1750882       | Total Running Time : | 2955 Hour 34 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 52                | %       | Manual ;           | 25           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 39.42             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 37.90             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 30.93             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 30.32             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 29.71             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 22.74             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 21.22             | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 1      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.70       | kN  | Min : | 3.28       | kN  | Punch NO.; | 71 |
| SD :                 | 6.00       | %   | Max : | 4.66       | kN  | Punch NO.; | 57 |
| Mean Main Pressure : | 28.59      | kN  | Min : | 25.25      | kN  | Punch NO.; | 77 |
| SD :                 | 5.20       | %   | Max : | 34.21      | kN  | Punch NO.; | 57 |
| Filling Depth :      | Increase ; | 174 | Times | Decrease ; | 170 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10                | RPM          |
| Pre Upper Punch Position ;    | 2.70    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.30              | mm           |
| Main Upper Punch Position ;   | 2.75    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.85              | mm           |
| Filling Depth :               | 4.06    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.802   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 18      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.48 | kN |             |    |
|                  | Max ;     | 0.92 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1750882            | Tablets |       |             |         |
| A.W.C Counter : | 1719061            | Tablets |       |             |         |
| HEP :           | 249 (0.0%)         | Tablets | LEP : | 1656 (0.1%) | Tablets |
| GOOD :          | 1717156<br>(99.9%) | Tablets |       |             |         |