

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | ms88737       | Batch NO. :          | KQM2510140              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 15 Min 58 Sec   |
| Total Counter :  | 1150768       | Total Running Time : | 2967 Hour 12 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 62                | %       | Manual ;           | 23           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 39.28             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 37.77             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 30.82             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 30.21             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 29.61             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.66             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 21.15             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.57  | kN    | Min :      | 3.20  | kN    | Punch NO.; | 19 |
|                      |            | 3.90  | %     | Max :      | 3.91  | kN    | Punch NO.; | 74 |
| Mean Main Pressure : | SD :       | 29.94 | kN    | Min :      | 26.17 | kN    | Punch NO.; | 9  |
|                      |            | 3.80  | %     | Max :      | 32.00 | kN    | Punch NO.; | 74 |
| Filling Depth :      | Increase ; | 56    | Times | Decrease ; | 65    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.69    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.31     | mm           |
| Main Upper Punch Position ;   | 2.69    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.91     | mm           |
| Filling Depth :               | 3.46    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.758   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.46 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1150768             | Tablets |       |            |         |
| A.W.C Counter : | 1145171             | Tablets |       |            |         |
| HEP :           | 6 (0.0%)            | Tablets | LEP : | 313 (0.0%) | Tablets |
| GOOD :          | 1144852<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | ms88737       | Batch NO. :          | KQM2510140              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 15 Min 58 Sec   |
| Total Counter :  | 1150768       | Total Running Time : | 2967 Hour 12 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 59                | %       | Manual ;           | 24           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 38.00             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 36.54             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 29.81             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 29.23             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 28.64             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 21.92             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 20.46             | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 1      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.59  | kN    | Min :      | 3.17  | kN    | Punch NO.; | 77 |
|                      |            | 5.10  | %     | Max :      | 4.01  | kN    | Punch NO.; | 48 |
| Mean Main Pressure : | SD :       | 29.19 | kN    | Min :      | 26.13 | kN    | Punch NO.; | 37 |
|                      |            | 4.40  | %     | Max :      | 32.45 | kN    | Punch NO.; | 48 |
| Filling Depth :      | Increase ; | 88    | Times | Decrease ; | 88    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 184800   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.70    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.00    | mm  | Pre Thickness ;  | 1.30     | mm           |
| Main Upper Punch Position ;   | 2.75    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.85     | mm           |
| Filling Depth :               | 3.39    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.758   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1150768             | Tablets |       |            |         |
| A.W.C Counter : | 1145171             | Tablets |       |            |         |
| HEP :           | 6 (0.0%)            | Tablets | LEP : | 313 (0.0%) | Tablets |
| GOOD :          | 1144852<br>(100.0%) | Tablets |       |            |         |