

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	06 Hour 31 Min 43 Sec
Total Counter :	1200584	Total Running Time :	2996 Hour 59 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	62	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	39.81	kN	No		
HEP :	25.0 %	38.28	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.24	kN	Yes 15	Times	
Ref :		30.63	kN			
LCP :	2.0 %	30.01	kN	Yes 15	Times	
LEP :	25.0 %	22.97	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.44	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.03	kN	Min :	3.59	kN	Punch NO.;	72
SD :	5.80	%	Max :	4.54	kN	Punch NO.;	43
Mean Main Pressure :	30.98	kN	Min :	27.41	kN	Punch NO.;	71
SD :	4.90	%	Max :	34.15	kN	Punch NO.;	4
Filling Depth :	Increase ;	432	Times	Decrease ;	398	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	3.90	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.50	mm	Main Thickness ;	0.80	mm
Filling Depth :	4.33	mm	Current Cam :	Cam B	
Main Air pressure :	0.648	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	48	Min			
Lower Punch Remain Time (S2):	20	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.28	kN		
	Max ;	0.38	kN	Max Punch ;	35
Lower Punch :	Average ;	0.42	kN		
	Max ;	1.51	kN	Max Punch ;	70
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1200584	Tablets			
A.W.C Counter :	1168611	Tablets			
HEP :	145 (0.0%)	Tablets	LEP :	363 (0.0%)	Tablets
GOOD :	1168103 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	06 Hour 31 Min 43 Sec
Total Counter :	1200584	Total Running Time :	2996 Hour 59 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	64	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	37.92	kN	No		
HEP :	25.0 %	36.46	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	29.75	kN	Yes 15	Times	
Ref :		29.17	kN			
LCP :	2.0 %	28.59	kN	Yes 15	Times	
LEP :	25.0 %	21.88	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	20.42	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.04	kN	Min :	3.68	kN	Punch NO.;	23
SD :	4.20	%	Max :	4.46	kN	Punch NO.;	56
Mean Main Pressure :	29.33	kN	Min :	27.10	kN	Punch NO.;	71
SD :	3.20	%	Max :	31.89	kN	Punch NO.;	19
Filling Depth :	Increase ;	119	Times	Decrease ;	103	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	22	RPM	Capacity :	203280	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.70	mm			
Pre Lower Punch Position ;	3.90	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	2.70	mm			
Main Lower Punch Position ;	3.50	mm	Main Thickness ;	0.80	mm
Filling Depth :	3.45	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.648	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	48	Min			
Lower Punch Remain Time (S2):	20	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.28	kN		
	Max ;	0.38	kN	Max Punch ;	35
Lower Punch :	Average ;	0.42	kN		
	Max ;	1.58	kN	Max Punch ;	70
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1200584	Tablets			
A.W.C Counter :	1168611	Tablets			
HEP :	145 (0.0%)	Tablets	LEP :	363 (0.0%)	Tablets
GOOD :	1168103 (100.0%)	Tablets			