

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	08 Hour 29 Min 04 Sec
Total Counter :	1600782	Total Running Time :	2998 Hour 57 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	74	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.48	kN	No		
HEP :	25.0 %	38.93	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.76	kN	Yes 15	Times	
Ref :		31.14	kN			
LCP :	2.0 %	30.52	kN	Yes 15	Times	
LEP :	25.0 %	23.36	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.80	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	4.18	kN	Min :	3.62	kN	Punch NO.;	5
SD :	6.40	%	Max :	4.98	kN	Punch NO.;	42
Mean Main Pressure :	31.34	kN	Min :	27.32	kN	Punch NO.;	6
SD :	5.40	%	Max :	35.44	kN	Punch NO.;	42
Filling Depth :	Increase ;	542	Times	Decrease ;	489	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	18	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.21	mm	Pre Thickness ;	1.21	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.84	mm	Main Thickness ;	0.84	mm
Filling Depth :	5.07	mm	Current Cam :	Cam B	
Main Air pressure :	0.631	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	22	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.40	kN		
	Max ;	0.51	kN	Max Punch ;	54
Lower Punch :	Average ;	0.42	kN		
	Max ;	0.99	kN	Max Punch ;	49
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1600782	Tablets			
A.W.C Counter :	1557693	Tablets			
HEP :	334 (0.0%)	Tablets	LEP :	1023 (0.1%)	Tablets
GOOD :	1556336 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510440
Print Interval :	50000 Tablets	Running Time :	08 Hour 29 Min 04 Sec
Total Counter :	1600782	Total Running Time :	2998 Hour 57 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	61	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	34.08	kN	No		
HEP :	25.0 %	32.77	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	26.74	kN	Yes 15	Times	
Ref :		26.22	kN			
LCP :	2.0 %	25.69	kN	Yes 15	Times	
LEP :	25.0 %	19.66	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	18.35	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.88	kN	Min :	3.35	kN	Punch NO.;	49
SD :	7.00	%	Max :	4.63	kN	Punch NO.;	22
Mean Main Pressure :	26.72	kN	Min :	23.19	kN	Punch NO.;	49
SD :	6.40	%	Max :	30.97	kN	Punch NO.;	19
Filling Depth :	Increase ;	262	Times	Decrease ;	236	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	221760	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.56	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.631	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	22	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.40	kN		
	Max ;	0.51	kN	Max Punch ;	54
Lower Punch :	Average ;	0.42	kN		
	Max ;	1.03	kN	Max Punch ;	49
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1600782	Tablets			
A.W.C Counter :	1557693	Tablets			
HEP :	334 (0.0%)	Tablets	LEP :	1023 (0.1%)	Tablets
GOOD :	1556336 (99.9%)	Tablets			