

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2510440              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 09 Min 54 Sec   |
| Total Counter :  | 1750876       | Total Running Time : | 2999 Hour 38 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 74                | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 10.0         | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 40.91             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 39.34             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 32.10             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 31.47             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 30.84             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 23.60             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 22.03             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 4.04       | kN  | Min : | 3.63       | kN  | Punch NO.; | 62 |
| SD :                 | 5.30       | %   | Max : | 4.50       | kN  | Punch NO.; | 53 |
| Mean Main Pressure : | 31.31      | kN  | Min : | 28.52      | kN  | Punch NO.; | 55 |
| SD :                 | 4.20       | %   | Max : | 34.16      | kN  | Punch NO.; | 4  |
| Filling Depth :      | Increase ; | 554 | Times | Decrease ; | 521 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 24      | RPM | Capacity :       | 221760   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.21    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.84    | mm  | Main Thickness ; | 0.84     | mm           |
| Filling Depth :               | 4.61    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.632   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 31      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.62 | kN |             |    |
|                  | Max ;     | 0.79 | kN | Max Punch ; | 77 |
| Lower Punch :    | Average ; | 0.49 | kN |             |    |
|                  | Max ;     | 1.14 | kN | Max Punch ; | 13 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1750876            | Tablets |       |             |         |
| A.W.C Counter : | 1707766            | Tablets |       |             |         |
| HEP :           | 335 (0.0%)         | Tablets | LEP : | 1028 (0.1%) | Tablets |
| GOOD :          | 1706403<br>(99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2510440              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 09 Min 54 Sec   |
| Total Counter :  | 1750876       | Total Running Time : | 2999 Hour 38 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 92                | %       | Manual ;           | 45    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 39.37             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 37.85             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 30.89             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 30.28             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 29.68             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 22.71             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 21.20             | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 1      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 4.32       | kN  | Min : | 3.80       | kN  | Punch NO.; | 77 |
| SD :                 | 6.10       | %   | Max : | 5.06       | kN  | Punch NO.; | 34 |
| Mean Main Pressure : | 29.88      | kN  | Min : | 26.77      | kN  | Punch NO.; | 77 |
| SD :                 | 4.90       | %   | Max : | 33.31      | kN  | Punch NO.; | 29 |
| Filling Depth :      | Increase ; | 366 | Times | Decrease ; | 341 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 24      | RPM | Capacity :       | 221760            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.20              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.85    | mm  | Main Thickness ; | 0.85              | mm           |
| Filling Depth :               | 3.90    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.632   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 31      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.62 | kN |             |    |
|                  | Max ;     | 0.80 | kN | Max Punch ; | 61 |
| Lower Punch :    | Average ; | 0.49 | kN |             |    |
|                  | Max ;     | 1.14 | kN | Max Punch ; | 13 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1750876            | Tablets |       |             |         |
| A.W.C Counter : | 1707766            | Tablets |       |             |         |
| HEP :           | 335 (0.0%)         | Tablets | LEP : | 1028 (0.1%) | Tablets |
| GOOD :          | 1706403<br>(99.9%) | Tablets |       |             |         |