

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510540
Print Interval :	50000 Tablets	Running Time :	08 Hour 11 Min 53 Sec
Total Counter :	1751038	Total Running Time :	3011 Hour 58 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	57	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	40.66	kN	No		
HEP :	25.0 %	39.09	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	31.90	kN	Yes 15	Times	
Ref :		31.28	kN			
LCP :	2.0 %	30.65	kN	Yes 15	Times	
LEP :	25.0 %	23.46	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	21.89	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.72	kN	Min :	3.28	kN	Punch NO.;	26
SD :	5.80	%	Max :	4.28	kN	Punch NO.;	38
Mean Main Pressure :	29.23	kN	Min :	25.82	kN	Punch NO.;	26
SD :	5.80	%	Max :	33.40	kN	Punch NO.;	38
Filling Depth :	Increase ;	1710	Times	Decrease ;	1673	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.73	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.747	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	24	Min			
Lower Punch Remain Time (S2):	3	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.46	kN		
	Max ;	0.59	kN	Max Punch ;	19
Lower Punch :	Average ;	0.29	kN		
	Max ;	0.83	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1751038	Tablets			
A.W.C Counter :	1725317	Tablets			
HEP :	559 (0.0%)	Tablets	LEP :	2115 (0.1%)	Tablets
GOOD :	1722643 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 40 MG,STKQ1B00
User ID :	pb97604	Batch NO. :	KQM2510540
Print Interval :	50000 Tablets	Running Time :	08 Hour 11 Min 53 Sec
Total Counter :	1751038	Total Running Time :	3011 Hour 58 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	62	%	Manual ;	45	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2750000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	36.59	kN	No		
HEP :	25.0 %	35.18	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	28.71	kN	Yes 15	Times	
Ref :		28.14	kN			
LCP :	2.0 %	27.58	kN	Yes 15	Times	
LEP :	25.0 %	21.11	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	19.70	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	1	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.93	kN	Min :	3.53	kN	Punch NO.;	70
SD :	4.80	%	Max :	4.47	kN	Punch NO.;	46
Mean Main Pressure :	27.98	kN	Min :	24.71	kN	Punch NO.;	70
SD :	4.60	%	Max :	30.72	kN	Punch NO.;	46
Filling Depth :	Increase ;	1521	Times	Decrease ;	1412	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	24	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.20	mm	Pre Thickness ;	1.20	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.85	mm	Main Thickness ;	0.85	mm
Filling Depth :	3.87	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.747	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	24	Min			
Lower Punch Remain Time (S2):	3	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.46	kN		
	Max ;	0.59	kN	Max Punch ;	35
Lower Punch :	Average ;	0.29	kN		
	Max ;	0.83	kN	Max Punch ;	40
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	1751038	Tablets			
A.W.C Counter :	1725394	Tablets			
HEP :	559 (0.0%)	Tablets	LEP :	2115 (0.1%)	Tablets
GOOD :	1722720 (99.8%)	Tablets			