

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2510740              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 21 Min 23 Sec   |
| Total Counter :  | 300150        | Total Running Time : | 3029 Hour 57 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 97                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 43.94             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 42.25             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 34.47             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 33.80             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 33.12             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 25.35             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 23.66             | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.23       | kN  | Min : | 2.17       | kN  | Punch NO.; | 13 |
| SD :                 | 1.60       | %   | Max : | 2.30       | kN  | Punch NO.; | 17 |
| Mean Main Pressure : | 34.00      | kN  | Min : | 31.20      | kN  | Punch NO.; | 3  |
| SD :                 | 3.90       | %   | Max : | 36.90      | kN  | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 394 | Times | Decrease ; | 300 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 24      | RPM | Capacity :       | 221760   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 23       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.85    | mm  | Main Thickness ; | 0.85     | mm           |
| Filling Depth :               | 3.84    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.746   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 36      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 30 |
| Lower Punch :    | Average ; | 0.28 | kN |             |    |
|                  | Max ;     | 0.70 | kN | Max Punch ; | 61 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 300150             | Tablets |       |           |         |
| A.W.C Counter : | 292007             | Tablets |       |           |         |
| HEP :           | 3 (0.0%)           | Tablets | LEP : | 25 (0.0%) | Tablets |
| GOOD :          | 291979<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40<br>MG,STKQ1B00 |
| User ID :        | rp98927       | Batch NO. :          | KQM2510740              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 21 Min 23 Sec   |
| Total Counter :  | 300150        | Total Running Time : | 3029 Hour 57 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 97                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 38.95             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 37.46             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 30.56             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 29.96             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 29.37             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 22.47             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 20.98             | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 1      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.26       | kN  | Min : | 2.19       | kN  | Punch NO.; | 18 |
| SD :                 | 1.40       | %   | Max : | 2.33       | kN  | Punch NO.; | 8  |
| Mean Main Pressure : | 29.62      | kN  | Min : | 25.78      | kN  | Punch NO.; | 27 |
| SD :                 | 4.30       | %   | Max : | 33.09      | kN  | Punch NO.; | 36 |
| Filling Depth :      | Increase ; | 393 | Times | Decrease ; | 365 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 24      | RPM | Capacity :       | 221760   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 23       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.20    | mm  | Pre Thickness ;  | 1.21     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.85    | mm  | Main Thickness ; | 0.86     | mm           |
| Filling Depth :               | 4.17    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.746   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 36      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 71 |
| Lower Punch :    | Average ; | 0.28 | kN |             |    |
|                  | Max ;     | 0.70 | kN | Max Punch ; | 61 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 300150             | Tablets |       |           |         |
| A.W.C Counter : | 292007             | Tablets |       |           |         |
| HEP :           | 3 (0.0%)           | Tablets | LEP : | 25 (0.0%) | Tablets |
| GOOD :          | 291979<br>(100.0%) | Tablets |       |           |         |