

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STHA1A00 |
| User ID :        | sj97904       | Batch NO. :          | HAM2501650              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 15 Min 11 Sec   |
| Total Counter :  | 500108        | Total Running Time : | 3081 Hour 48 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 98       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 750000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 50       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 %      | 11.89    | kN      | No                 |              |           |
| HEP :                              | 15.0 %      | 11.12    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 9.81     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.67     | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 9.52     | kN      | Yes 15             | Times        |           |
| LEP :                              | 15.0 %      | 8.22     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 %      | 7.44     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 2.05       | kN   | Min : | 2.04       | kN   | Punch NO.; | 4  |
| SD :                 | 0.30       | %    | Max : | 2.08       | kN   | Punch NO.; | 7  |
| Mean Main Pressure : | 9.62       | kN   | Min : | 9.14       | kN   | Punch NO.; | 52 |
| SD :                 | 2.30       | %    | Max : | 10.15      | kN   | Punch NO.; | 17 |
| Filling Depth :      | Increase ; | 2131 | Times | Decrease ; | 2003 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.70    | mm  | Main Thickness ; | 2.70     | mm           |
| Filling Depth :               | 7.29    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.762   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 6  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500108            | Tablets |       |            |         |
| A.W.C Counter : | 497238            | Tablets |       |            |         |
| HEP :           | 615 (0.1%)        | Tablets | LEP : | 233 (0.0%) | Tablets |
| GOOD :          | 496390<br>(99.8%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STHA1A00 |
| User ID :        | sj97904       | Batch NO. :          | HAM2501650              |
| Print Interval : | 50000 Tablets | Running Time :       | 05 Hour 15 Min 11 Sec   |
| Total Counter :  | 500108        | Total Running Time : | 3081 Hour 48 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 98                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 750000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 50                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 23.0 % | 12.46             | kN      | No                 |              |           |
| HEP :                              | 15.0 % | 11.65             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 10.29             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 10.13             | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 9.98              | kN      | Yes 15             | Times        |           |
| LEP :                              | 15.0 % | 8.61              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 23.0 % | 7.80              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 2.08       | kN   | Min : | 2.06       | kN   | Punch NO.; | 42 |
| SD :                 | 0.60       | %    | Max : | 2.11       | kN   | Punch NO.; | 2  |
| Mean Main Pressure : | 10.29      | kN   | Min : | 9.79       | kN   | Punch NO.; | 4  |
| SD :                 | 2.50       | %    | Max : | 10.99      | kN   | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 2499 | Times | Decrease ; | 2504 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.85    | mm  | Main Thickness ; | 2.85     | mm           |
| Filling Depth :               | 7.35    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.762   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 6  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500108            | Tablets |       |            |         |
| A.W.C Counter : | 497238            | Tablets |       |            |         |
| HEP :           | 615 (0.1%)        | Tablets | LEP : | 233 (0.0%) | Tablets |
| GOOD :          | 496390<br>(99.8%) | Tablets |       |            |         |