

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRALINE 50mg<br>TABS |
| User ID :        | rp98927       | Batch NO. :          | SED-TAB-XV-25-001A      |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 56 Min 47 Sec   |
| Total Counter :  | 450350        | Total Running Time : | 3187 Hour 33 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 16       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 700000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 3.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.90     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 11.26    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 10.83    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %       | 8.79     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 8.66     | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 8.53     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 6.50     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 6.07     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.44       | kN  | Min : | 1.40       | kN  | Punch NO.; | 58 |
| SD :                 | 0.90       | %   | Max : | 1.46       | kN  | Punch NO.; | 2  |
| Mean Main Pressure : | 8.68       | kN  | Min : | 8.28       | kN  | Punch NO.; | 6  |
| SD :                 | 1.90       | %   | Max : | 9.19       | kN  | Punch NO.; | 62 |
| Filling Depth :      | Increase ; | 215 | Times | Decrease ; | 263 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.70    | mm  | Main Thickness ; | 2.71     | mm           |
| Filling Depth :               | 5.90    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.776   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 65 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 450350             | Tablets |       |           |         |
| A.W.C Counter : | 439943             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 42 (0.0%) | Tablets |
| GOOD :          | 439901<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRALINE 50mg<br>TABS |
| User ID :        | rp98927       | Batch NO. :          | SED-TAB-XV-25-001A      |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 56 Min 47 Sec   |
| Total Counter :  | 450350        | Total Running Time : | 3187 Hour 33 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 17                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 700000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 30                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 3.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.90              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 10.70             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.29             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 1.5 %  | 8.35              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.23              | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 8.11              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.17              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 5.76              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.29       | kN  | Min : | 1.25       | kN  | Punch NO.; | 35 |
| SD :                 | 1.40       | %   | Max : | 1.33       | kN  | Punch NO.; | 23 |
| Mean Main Pressure : | 8.17       | kN  | Min : | 7.69       | kN  | Punch NO.; | 69 |
| SD :                 | 3.10       | %   | Max : | 9.34       | kN  | Punch NO.; | 62 |
| Filling Depth :      | Increase ; | 112 | Times | Decrease ; | 139 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.90    | mm  | Main Thickness ; | 2.91     | mm           |
| Filling Depth :               | 5.97    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.776   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 4  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 35 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 450350             | Tablets |       |           |         |
| A.W.C Counter : | 439943             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 42 (0.0%) | Tablets |
| GOOD :          | 439901<br>(100.0%) | Tablets |       |           |         |