

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG,STIU1B00   |
| User ID :        | cb83530       | Batch NO. :          | IUM2501240            |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 28 Min 36 Sec |
| Total Counter :  | 1650576       | Total Running Time : | 3205 Hour 48 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 180      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2025000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20       | Min     | Supply Time ;      | 1.1          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.8          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 26.18    | kN      | No                 |              |           |
| HEP :                              | 24.0 %      | 24.97    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 20.54    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 20.14    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.73    | kN      | Yes 20             | Times        |           |
| LEP :                              | 24.0 %      | 15.30    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 14.10    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.57       | kN  | Min : | 2.48       | kN  | Punch NO.; | 2  |
| SD :                 | 1.70       | %   | Max : | 2.69       | kN  | Punch NO.; | 25 |
| Mean Main Pressure : | 20.17      | kN  | Min : | 18.64      | kN  | Punch NO.; | 41 |
| SD :                 | 4.30       | %   | Max : | 23.52      | kN  | Punch NO.; | 25 |
| Filling Depth :      | Increase ; | 253 | Times | Decrease ; | 237 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 32       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.75    | mm  | Pre Thickness ;  | 2.76     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.34    | mm  | Main Thickness ; | 2.35     | mm           |
| Filling Depth :               | 7.05    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.760   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 1650576            | Tablets |       |           |         |
| A.W.C Counter : | 1644682            | Tablets |       |           |         |
| HEP :           | 1312 (0.1%)        | Tablets | LEP : | 64 (0.0%) | Tablets |
| GOOD :          | 1643306<br>(99.9%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBU 400 MG,STIU1B00   |
| User ID :        | cb83530       | Batch NO. :          | IUM2501240            |
| Print Interval : | 50000 Tablets | Running Time :       | 14 Hour 28 Min 36 Sec |
| Total Counter :  | 1650576       | Total Running Time : | 3205 Hour 48 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 180      | %       | Manual ;           | 21           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2025000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 20       | Min     | Supply Time ;      | 1.1          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.8          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 23.16    | kN      | No                 |              |           |
| HEP :                              | 24.0 %      | 22.09    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 18.17    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 17.81    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 17.46    | kN      | Yes 20             | Times        |           |
| LEP :                              | 24.0 %      | 13.54    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 12.47    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.40           | kN    | Min :      | 2.31  | kN    | Punch NO.; | 2  |
| SD :                 | 2.00           | %     | Max :      | 2.52  | kN    | Punch NO.; | 25 |
| Mean Main Pressure : | 17.77          | kN    | Min :      | 16.13 | kN    | Punch NO.; | 17 |
| SD :                 | 4.40           | %     | Max :      | 20.62 | kN    | Punch NO.; | 25 |
| Filling Depth :      | Increase ; 193 | Times | Decrease ; | 175   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 32       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.80    | mm  | Pre Thickness ;  | 2.81     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.45    | mm  | Main Thickness ; | 2.46     | mm           |
| Filling Depth :               | 6.99    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.760   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 44 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                 |         |       |           |         |
|-----------------|-----------------|---------|-------|-----------|---------|
| Total Counter : | 1650576         | Tablets |       |           |         |
| A.W.C Counter : | 1644682         | Tablets |       |           |         |
| HEP :           | 1312 (0.1%)     | Tablets | LEP : | 64 (0.0%) | Tablets |
| GOOD :          | 1643306 (99.9%) | Tablets |       |           |         |