

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	ATORV 80 MG,STJB1D00
User ID :	vg96371	Batch NO. :	JBM2502580
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 00 Sec
Total Counter :	50008	Total Running Time :	3276 Hour 40 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	319	%	Manual ;	160	RPM
Filling Cam :		Cam E 16.5- 22mm		Cam Range :	Min :	14.50 mm
Target Quantity :		750000	Tablets		Max :	20.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	4.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	18.31	kN	No		
HEP :	25.0 %	17.60	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	14.36	kN	Yes	15	Times
Ref :		14.08	kN			
LCP :	2.0 %	13.80	kN	Yes	15	Times
LEP :	25.0 %	10.56	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	9.86	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.81	kN	Min :	2.64	kN	Punch NO.;	39
SD :	3.30	%	Max :	3.06	kN	Punch NO.;	48
Mean Main Pressure :	14.12	kN	Min :	13.07	kN	Punch NO.;	1
SD :	4.50	%	Max :	16.15	kN	Punch NO.;	48
Filling Depth :	Increase ;	296	Times	Decrease ;	236	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	48	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	7.40	mm	Main Thickness ;	4.40	mm
Filling Depth :	16.19	mm	Current Cam :	Cam E 16.5- 22mm	
Main Air pressure :	0.749	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	35	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.01	kN		
	Max ;	0.05	kN	Max Punch ;	35
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.21	kN	Max Punch ;	25
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50008	Tablets			
A.W.C Counter :	46404	Tablets			
HEP :	112 (0.2%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	46292 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	ATORV 80 MG,STJB1D00
User ID :	vg96371	Batch NO. :	JBM2502580
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 00 Sec
Total Counter :	50008	Total Running Time :	3276 Hour 40 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	335	%	Manual ;	150	RPM
Filling Cam :		Cam E 16.5- 22mm		Cam Range :	Min :	14.50 mm
Target Quantity :		750000	Tablets		Max :	20.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		50	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	4.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		3.00	kN			
Ejecting Force High Limit :		4.90	kN			
				Stop Condition ;		
HSP :	30.0 %	17.50	kN	No		
HEP :	25.0 %	16.83	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	13.73	kN	Yes	15	Times
Ref :		13.46	kN			
LCP :	2.0 %	13.20	kN	Yes	15	Times
LEP :	25.0 %	10.10	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	9.43	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.83	kN	Min :	2.72	kN	Punch NO.;	2
SD :	2.20	%	Max :	3.06	kN	Punch NO.;	22
Mean Main Pressure :	13.50	kN	Min :	12.65	kN	Punch NO.;	38
SD :	2.90	%	Max :	14.92	kN	Punch NO.;	22
Filling Depth :	Increase ;	146	Times	Decrease ;	166	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	8.99	mm	Pre Thickness ;	6.00	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	7.55	mm	Main Thickness ;	4.55	mm
Filling Depth :	15.64	mm	Current Cam :	Cam E 16.5- 22mm	
Main Air pressure :	0.749	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	35	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.01	kN		
	Max ;	0.04	kN	Max Punch ;	35
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.21	kN	Max Punch ;	25
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	50008	Tablets			
A.W.C Counter :	46457	Tablets			
HEP :	112 (0.2%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	46345 (99.8%)	Tablets			