

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ATORV 80<br>MG,STJB1D00 |
| User ID :        | pb97604       | Batch NO. :          | JBM2502680              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 23 Min 59 Sec   |
| Total Counter :  | 250062        | Total Running Time : | 3285 Hour 48 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                        |         |                    |       |                      |
|------------------------------------|--------|------------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 380                    | %       | Manual ;           | 160   | RPM                  |
| Filling Cam :                      |        | Cam E<br>16.5-<br>22mm |         | Cam Range :        | Min : | 14.50 mm             |
| Target Quantity :                  |        | 750000                 | Tablets |                    | Max : | 20.00 mm             |
| Air Pressure Low Limit :           |        | 0.400                  | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                        |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                     | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                     | Min     | Supply Time ;      | 4.0   | Sec                  |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                      | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00                   | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00                   | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00                   | kN      |                    |       |                      |
|                                    |        |                        |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 16.65                  | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 16.01                  | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 13.06                  | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 12.81                  | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 12.55                  | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 9.60                   | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 8.96                   | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                        |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                        |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.99       | kN  | Min : | 2.78       | kN  | Punch NO.; | 2  |
| SD :                 | 3.70       | %   | Max : | 3.29       | kN  | Punch NO.; | 3  |
| Mean Main Pressure : | 12.97      | kN  | Min : | 11.82      | kN  | Punch NO.; | 38 |
| SD :                 | 4.80       | %   | Max : | 14.61      | kN  | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 447 | Times | Decrease ; | 403 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 65                     | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.80    | mm  | Pre Thickness ;  | 5.80                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.40    | mm  | Main Thickness ; | 4.40                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 15.52   | mm  | Current Cam :    | Cam E<br>16.5-<br>22mm |              |
| Main Air pressure :           | 0.743   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 6       | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 25 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 250062             | Tablets |       |        |         |
| A.W.C Counter : | 243339             | Tablets |       |        |         |
| HEP :           | 35 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 243304<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ATORV 80<br>MG,STJB1D00 |
| User ID :        | pb97604       | Batch NO. :          | JBM2502680              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 23 Min 59 Sec   |
| Total Counter :  | 250062        | Total Running Time : | 3285 Hour 48 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                        |         |                    |       |                      |
|------------------------------------|--------|------------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 380                    | %       | Manual ;           | 150   | RPM                  |
| Filling Cam :                      |        | Cam E<br>16.5-<br>22mm |         | Cam Range :        | Min : | 14.50 mm             |
| Target Quantity :                  |        | 750000                 | Tablets |                    | Max : | 20.00 mm             |
| Air Pressure Low Limit :           |        | 0.400                  | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                        |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 50                     | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 30                     | Min     | Supply Time ;      | 4.0   | Sec                  |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                      | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00                   | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00                   | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00                   | kN      |                    |       |                      |
|                                    |        |                        |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 17.43                  | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 16.76                  | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 13.68                  | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 13.41                  | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 13.14                  | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.06                  | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.39                   | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                        |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                        |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.03       | kN  | Min : | 2.88       | kN  | Punch NO.; | 39 |
| SD :                 | 2.50       | %   | Max : | 3.25       | kN  | Punch NO.; | 48 |
| Mean Main Pressure : | 13.44      | kN  | Min : | 12.77      | kN  | Punch NO.; | 15 |
| SD :                 | 3.10       | %   | Max : | 14.77      | kN  | Punch NO.; | 48 |
| Filling Depth :      | Increase ; | 184 | Times | Decrease ; | 220 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 65                     | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.80    | mm  | Pre Thickness ;  | 5.81                   | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.55    | mm  | Main Thickness ; | 4.55                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 15.45   | mm  | Current Cam :    | Cam E<br>16.5-<br>22mm |              |
| Main Air pressure :           | 0.743   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 6       | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 16      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.15 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 25 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 250062             | Tablets |       |        |         |
| A.W.C Counter : | 243339             | Tablets |       |        |         |
| HEP :           | 35 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 243304<br>(100.0%) | Tablets |       |        |         |