

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 750MG<br>STHA1B00 |
| User ID :        | mk98268       | Batch NO. :          | HAM2401075              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 27 Min 29 Sec   |
| Total Counter :  | 327944        | Total Running Time : | 322 Hour 14 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 255               | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 500000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 11.27             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.84             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.84              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.67              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.50              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.50              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.07              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None   |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 1.88       | kN   | Min : | 1.81       | kN   | Punch NO.; | 9  |
| SD :                 | 2.00       | %    | Max : | 1.95       | kN   | Punch NO.; | 32 |
| Mean Main Pressure : | 8.63       | kN   | Min : | 8.23       | kN   | Punch NO.; | 23 |
| SD :                 | 2.50       | %    | Max : | 9.11       | kN   | Punch NO.; | 13 |
| Filling Depth :      | Increase ; | 1585 | Times | Decrease ; | 1350 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 38                | RPM          |
| Pre Upper Punch Position ;    | 2.49    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.80    | mm  | Pre Thickness ;  | 5.31              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.55    | mm  | Main Thickness ; | 3.05              | mm           |
| Filling Depth :               | 9.56    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.820   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 51      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 33 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 15 |

### Tablet Data (N/A)

|                 |                    |         |       |          |         |
|-----------------|--------------------|---------|-------|----------|---------|
| Total Counter : | 327944             | Tablets |       |          |         |
| A.W.C Counter : | 326427             | Tablets |       |          |         |
| HEP :           | 3 (0.0%)           | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 326423<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 750MG<br>STHA1B00 |
| User ID :        | mk98268       | Batch NO. :          | HAM2401075              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 27 Min 29 Sec   |
| Total Counter :  | 327944        | Total Running Time : | 322 Hour 14 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 255               | %       | Manual ;           | 59           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 500000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 12.64             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 12.16             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 9.92              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.73              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 9.53              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 7.30              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.81              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |      |       |            |      |            |    |
|----------------------|------------|------|-------|------------|------|------------|----|
| Mean Pre Pressure :  | 1.92       | kN   | Min : | 1.86       | kN   | Punch NO.; | 14 |
| SD :                 | 2.30       | %    | Max : | 2.05       | kN   | Punch NO.; | 50 |
| Mean Main Pressure : | 9.84       | kN   | Min : | 8.88       | kN   | Punch NO.; | 31 |
| SD :                 | 4.40       | %    | Max : | 11.56      | kN   | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 1420 | Times | Decrease ; | 1277 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 38                | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.80    | mm  | Pre Thickness ;  | 5.30              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.55    | mm  | Main Thickness ; | 3.05              | mm           |
| Filling Depth :               | 9.68    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.820   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 51      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 33 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 9  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 15 |

### Tablet Data (N/A)

|                 |                    |         |       |          |         |
|-----------------|--------------------|---------|-------|----------|---------|
| Total Counter : | 327944             | Tablets |       |          |         |
| A.W.C Counter : | 326427             | Tablets |       |          |         |
| HEP :           | 3 (0.0%)           | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 326423<br>(100.0%) | Tablets |       |          |         |