

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MGSTVI1A0<br>0 |
| User ID :        | vg96371       | Batch NO. :          | VIM2401820              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 33 Min 42 Sec   |
| Total Counter :  | 1156552       | Total Running Time : | 330 Hour 33 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 76                | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 1740000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40                | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 13.48             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 12.96             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 10.57             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 10.37             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.16             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 7.78              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 7.26              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 0.65         | kN    | Min :      | 0.63  | kN    | Punch NO.; | 29 |
| SD :                 | 1.20         | %     | Max :      | 0.67  | kN    | Punch NO.; | 60 |
| Mean Main Pressure : | 10.36        | kN    | Min :      | 9.71  | kN    | Punch NO.; | 19 |
| SD :                 | 2.40         | %     | Max :      | 11.06 | kN    | Punch NO.; | 50 |
| Filling Depth :      | Increase ; 2 | Times | Decrease ; | 21    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.19    | mm  | Pre Thickness ;  | 4.69     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 1.40     | mm           |
| Filling Depth :               | 5.67    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 52      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 65 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1156552            | Tablets |       |            |         |
| A.W.C Counter : | 1136057            | Tablets |       |            |         |
| HEP :           | 513 (0.0%)         | Tablets | LEP : | 562 (0.0%) | Tablets |
| GOOD :          | 1134982<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MGSTVI1A0<br>0 |
| User ID :        | vg96371       | Batch NO. :          | VIM2401820              |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 33 Min 42 Sec   |
| Total Counter :  | 1156552       | Total Running Time : | 330 Hour 33 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 75                | %       | Manual ;           | 59           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 1740000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 40                | Min     | Supply Time ;      | 0.2          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 14.05             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 13.51             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 11.03             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 10.81             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 10.59             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 8.11              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 7.57              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 0.76       | kN | Min : | 0.73       | kN | Punch NO.; | 14 |
| SD :                 | 1.70       | %  | Max : | 0.79       | kN | Punch NO.; | 52 |
| Mean Main Pressure : | 10.78      | kN | Min : | 10.26      | kN | Punch NO.; | 19 |
| SD :                 | 2.10       | %  | Max : | 11.34      | kN | Punch NO.; | 7  |
| Filling Depth :      | Increase ; | 13 | Times | Decrease ; | 19 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15                | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.39    | mm  | Pre Thickness ;  | 4.89              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.93    | mm  | Main Thickness ; | 1.43              | mm           |
| Filling Depth :               | 5.89    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.802   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 52      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 4       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 5  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 1156552            | Tablets |       |            |         |
| A.W.C Counter : | 1136057            | Tablets |       |            |         |
| HEP :           | 513 (0.0%)         | Tablets | LEP : | 562 (0.0%) | Tablets |
| GOOD :          | 1134982<br>(99.9%) | Tablets |       |            |         |