

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MG<br>STV11A00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500120              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 49 Min 24 Sec   |
| Total Counter :  | 700434        | Total Running Time : | 346 Hour 53 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 70                | %       | Manual ;           | 45    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 1740000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.3   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 60                | Min     | Supply Time ;      | 0.3   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 12.71             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 12.23             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 9.98              | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 9.78              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 9.59              | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 7.34              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 6.85              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 0.91       | kN | Min : | 0.88       | kN | Punch NO.; | 15 |
| SD :                 | 1.60       | %  | Max : | 0.94       | kN | Punch NO.; | 2  |
| Mean Main Pressure : | 9.88       | kN | Min : | 9.32       | kN | Punch NO.; | 9  |
| SD :                 | 2.20       | %  | Max : | 10.40      | kN | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 4  | Times | Decrease ; | 11 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.90    | mm  | Pre Thickness ;  | 4.40     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 1.40     | mm           |
| Filling Depth :               | 5.53    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.789   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 16      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 38      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 28 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 700434            | Tablets |       |            |         |
| A.W.C Counter : | 689754            | Tablets |       |            |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 464 (0.1%) | Tablets |
| GOOD :          | 689290<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO200MG<br>STV11A00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500120              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 49 Min 24 Sec   |
| Total Counter :  | 700434        | Total Running Time : | 346 Hour 53 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 76                | %       | Manual ;           | 59           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 1740000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60                | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 13.76             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 13.23             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 10.79             | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 10.58             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 10.37             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 7.94              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 7.41              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.03       | kN | Min : | 1.00       | kN | Punch NO.; | 11 |
| SD :                 | 1.40       | %  | Max : | 1.06       | kN | Punch NO.; | 16 |
| Mean Main Pressure : | 10.51      | kN | Min : | 9.70       | kN | Punch NO.; | 13 |
| SD :                 | 2.50       | %  | Max : | 10.99      | kN | Punch NO.; | 25 |
| Filling Depth :      | Increase ; | 5  | Times | Decrease ; | 17 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19                | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.10    | mm  | Pre Thickness ;  | 4.60              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.94    | mm  | Main Thickness ; | 1.44              | mm           |
| Filling Depth :               | 5.75    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.789   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 16      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 38      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 24 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 28 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 700434            | Tablets |       |            |         |
| A.W.C Counter : | 689819            | Tablets |       |            |         |
| HEP :           | 0 (0%)            | Tablets | LEP : | 464 (0.1%) | Tablets |
| GOOD :          | 689355<br>(99.9%) | Tablets |       |            |         |