

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | vg96371       | Batch NO. :          | VIM2500140              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 25 Min 56 Sec   |
| Total Counter :  | 50024         | Total Running Time : | 363 Hour 50 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 70       | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60       | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.06    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.33    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 14.96    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 14.66    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.37    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 11.00    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 10.26    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.58         | kN    | Min :        | 2.48  | kN | Punch NO.; | 24 |
| SD :                 | 1.50         | %     | Max :        | 2.65  | kN | Punch NO.; | 43 |
| Mean Main Pressure : | 14.74        | kN    | Min :        | 14.14 | kN | Punch NO.; | 9  |
| SD :                 | 1.80         | %     | Max :        | 15.44 | kN | Punch NO.; | 30 |
| Filling Depth :      | Increase ; 2 | Times | Decrease ; 5 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 117000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.44    | mm  | Main Thickness ; | 1.94     | mm           |
| Filling Depth :               | 7.26    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.794   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 65      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 50024            | Tablets |       |            |         |
| A.W.C Counter : | 45465            | Tablets |       |            |         |
| HEP :           | 43 (0.1%)        | Tablets | LEP : | 125 (0.3%) | Tablets |
| GOOD :          | 45297<br>(99.6%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | vg96371       | Batch NO. :          | VIM2500140              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 25 Min 56 Sec   |
| Total Counter :  | 50024         | Total Running Time : | 363 Hour 50 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 70                | %       | Manual ;           | 59           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 1500000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.3          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.3          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 19.17             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 18.44             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 15.04             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 14.75             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 14.45             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 11.06             | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 10.32             | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.37       | kN | Min : | 2.28       | kN | Punch NO.; | 32 |
| SD :                 | 1.90       | %  | Max : | 2.51       | kN | Punch NO.; | 46 |
| Mean Main Pressure : | 14.66      | kN | Min : | 13.99      | kN | Punch NO.; | 22 |
| SD :                 | 1.90       | %  | Max : | 15.42      | kN | Punch NO.; | 30 |
| Filling Depth :      | Increase ; | 1  | Times | Decrease ; | 0  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 117000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.52    | mm  | Main Thickness ; | 2.02     | mm           |
| Filling Depth :               | 7.56    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.794   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 65      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 50024            | Tablets |       |            |         |
| A.W.C Counter : | 45465            | Tablets |       |            |         |
| HEP :           | 43 (0.1%)        | Tablets | LEP : | 125 (0.3%) | Tablets |
| GOOD :          | 45297<br>(99.6%) | Tablets |       |            |         |