

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500140              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 18 Min 18 Sec   |
| Total Counter :  | 600380        | Total Running Time : | 366 Hour 42 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 89                | %       | Manual ;           | 45           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 1500000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.70             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.94             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 15.46             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 15.16             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.85             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 11.37             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 10.61             | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.46         | kN    | Min :      | 2.38  | kN    | Punch NO.; | 16 |
| SD :                 | 1.50         | %     | Max :      | 2.55  | kN    | Punch NO.; | 43 |
| Mean Main Pressure : | 15.25        | kN    | Min :      | 14.53 | kN    | Punch NO.; | 17 |
| SD :                 | 1.90         | %     | Max :      | 15.83 | kN    | Punch NO.; | 58 |
| Filling Depth :      | Increase ; 4 | Times | Decrease ; | 13    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.44    | mm  | Main Thickness ; | 1.94     | mm           |
| Filling Depth :               | 7.26    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.827   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 60      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 22      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 11 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 12 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 600380             | Tablets |       |            |         |
| A.W.C Counter : | 590866             | Tablets |       |            |         |
| HEP :           | 43 (0.0%)          | Tablets | LEP : | 125 (0.0%) | Tablets |
| GOOD :          | 590698<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG<br>STV11B00 |
| User ID :        | pb97604       | Batch NO. :          | VIM2500140              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 18 Min 18 Sec   |
| Total Counter :  | 600380        | Total Running Time : | 366 Hour 42 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 89                | %       | Manual ;           | 59           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 1500000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60                | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 20.02             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 19.25             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 15.71             | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 15.40             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 15.09             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 11.55             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 10.78             | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.39       | kN | Min : | 2.30       | kN | Punch NO.; | 63 |
| SD :                 | 2.00       | %  | Max : | 2.57       | kN | Punch NO.; | 46 |
| Mean Main Pressure : | 15.52      | kN | Min : | 14.81      | kN | Punch NO.; | 14 |
| SD :                 | 2.10       | %  | Max : | 16.54      | kN | Punch NO.; | 30 |
| Filling Depth :      | Increase ; | 1  | Times | Decrease ; | 12 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 195000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22                | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.50              | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.52    | mm  | Main Thickness ; | 2.02              | mm           |
| Filling Depth :               | 7.46    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.827   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 60      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 22      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 11 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 12 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 600380             | Tablets |       |            |         |
| A.W.C Counter : | 590866             | Tablets |       |            |         |
| HEP :           | 43 (0.0%)          | Tablets | LEP : | 125 (0.0%) | Tablets |
| GOOD :          | 590698<br>(100.0%) | Tablets |       |            |         |