

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	AMIOD 200MG STAI3A00
User ID :	rp98927	Batch NO. :	AIM2500120
Print Interval :	50000 Tablets	Running Time :	00 Hour 58 Min 27 Sec
Total Counter :	100044	Total Running Time :	479 Hour 7 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	100	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		1500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	24.06	kN	No		
HEP :	20.0 %	22.21	kN	Yes	10 Rotations	10 Tablets
HCP :	2.0 %	18.88	kN	Yes	15 Times	
Ref :		18.51	kN			
LCP :	2.0 %	18.14	kN	Yes	15 Times	
LEP :	20.0 %	14.81	kN	Yes	10 Rotations	10 Tablets
LSP :	30.0 %	12.96	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.36	kN	Min :	1.35	kN	Punch NO.;	11
SD :	0.60	%	Max :	1.38	kN	Punch NO.;	7
Mean Main Pressure :	18.06	kN	Min :	16.73	kN	Punch NO.;	3
SD :	3.80	%	Max :	19.74	kN	Punch NO.;	19
Filling Depth :	Increase ;	74	Times	Decrease ;	56	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	18	RPM	Capacity :	114480	Tablets/hour
Feeder : Status ;	Auto		Speed ;	18	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	7.10	mm	Pre Thickness ;	4.11	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	3.95	mm	Main Thickness ;	0.96	mm
Filling Depth :	5.95	mm	Current Cam :	Cam B	
Main Air pressure :	0.752	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	25	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.30	kN		
	Max ;	0.37	kN	Max Punch ;	53
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.16	kN	Max Punch ;	26
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	41

Tablet Data (N/A)

Total Counter :	100044	Tablets			
A.W.C Counter :	90493	Tablets			
HEP :	0 (0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	90492 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	AMIOD 200MG STAI3A00
User ID :	rp98927	Batch NO. :	AIM2500120
Print Interval :	50000 Tablets	Running Time :	00 Hour 58 Min 27 Sec
Total Counter :	100044	Total Running Time :	479 Hour 7 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	100	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		1500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	22.43	kN	No		
HEP :	20.0 %	20.70	kN	Yes	10	Rotations 10 Tablets
HCP :	2.0 %	17.60	kN	Yes	15	Times
Ref :		17.25	kN			
LCP :	2.0 %	16.91	kN	Yes	15	Times
LEP :	20.0 %	13.80	kN	Yes	10	Rotations 10 Tablets
LSP :	30.0 %	12.08	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.10	kN	Min :	1.07	kN	Punch NO.;	29
SD :	0.80	%	Max :	1.13	kN	Punch NO.;	8
Mean Main Pressure :	17.27	kN	Min :	16.19	kN	Punch NO.;	22
SD :	3.10	%	Max :	18.49	kN	Punch NO.;	19
Filling Depth :	Increase ;	42	Times	Decrease ;	24	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	18	RPM	Capacity :	114480	Tablets/hour
Feeder : Status ;	Auto		Speed ;	18	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	7.00	mm	Pre Thickness ;	4.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.01	mm	Main Thickness ;	1.02	mm
Filling Depth :	6.09	mm	Current Cam :	Cam B	
Main Air pressure :	0.752	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	25	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.29	kN		
	Max ;	0.38	kN	Max Punch ;	53
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.16	kN	Max Punch ;	26
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	41

Tablet Data (N/A)

Total Counter :	100044	Tablets			
A.W.C Counter :	90546	Tablets			
HEP :	0 (0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	90545 (100.0%)	Tablets			