

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO800mg<br>STV11C00 |
| User ID :        | yu98249       | Batch NO. :          | VIM2500380              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 55 Min 36 Sec   |
| Total Counter :  | 100040        | Total Running Time : | 505 Hour 2 Min          |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 149               | %       | Manual ;           | 35           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 750000            | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 24.00             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 23.07             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 18.83             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 18.46             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 18.09             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 13.84             | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 12.92             | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 20.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.30          | kN    | Min :      | 2.19  | kN    | Punch NO.; | 24 |
| SD :                 | 2.20          | %     | Max :      | 2.42  | kN    | Punch NO.; | 49 |
| Mean Main Pressure : | 18.54         | kN    | Min :      | 17.93 | kN    | Punch NO.; | 42 |
| SD :                 | 1.40          | %     | Max :      | 19.21 | kN    | Punch NO.; | 27 |
| Filling Depth :      | Increase ; 24 | Times | Decrease ; | 75    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 30       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.40    | mm  | Pre Thickness ;  | 4.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.77    | mm  | Main Thickness ; | 2.78     | mm           |
| Filling Depth :               | 10.03   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.752   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 39      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 9  |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 100040           | Tablets |       |            |         |
| A.W.C Counter : | 93301            | Tablets |       |            |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 102 (0.1%) | Tablets |
| GOOD :          | 93199<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO800mg<br>STV11C00 |
| User ID :        | yu98249       | Batch NO. :          | VIM2500380              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 55 Min 36 Sec   |
| Total Counter :  | 100040        | Total Running Time : | 505 Hour 2 Min          |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 149       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 750000    | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80        | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 23.85     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 22.94     | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 18.72     | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 18.35     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 17.98     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 13.76     | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 12.84     | kN      | No                 |              |           |
| SD limit :                         | 20.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 20.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.78       | kN | Min : | 2.69       | kN | Punch NO.; | 47 |
| SD :                 | 1.90       | %  | Max : | 2.91       | kN | Punch NO.; | 16 |
| Mean Main Pressure : | 18.49      | kN | Min : | 18.00      | kN | Punch NO.; | 18 |
| SD :                 | 1.80       | %  | Max : | 19.60      | kN | Punch NO.; | 27 |
| Filling Depth :      | Increase ; | 24 | Times | Decrease ; | 20 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 30       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.40    | mm  | Pre Thickness ;  | 4.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.83    | mm  | Main Thickness ; | 2.84     | mm           |
| Filling Depth :               | 9.98    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.752   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 39      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 9  |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 100040           | Tablets |       |            |         |
| A.W.C Counter : | 93301            | Tablets |       |            |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 102 (0.1%) | Tablets |
| GOOD :          | 93199<br>(99.9%) | Tablets |       |            |         |