

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | vg96371       | Batch NO. :          | VIM2500340            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 33 Min 56 Sec |
| Total Counter :  | 550304        | Total Running Time : | 514 Hour 0 Min        |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 75       | %       | Manual ;           | 35           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.87    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.14    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 14.80    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.51    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.22    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 10.88    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 10.16    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.14         | kN    | Min :      | 2.04  | kN    | Punch NO.; | 14 |
| SD :                 | 2.20         | %     | Max :      | 2.26  | kN    | Punch NO.; | 6  |
| Mean Main Pressure : | 14.42        | kN    | Min :      | 14.00 | kN    | Punch NO.; | 54 |
| SD :                 | 1.40         | %     | Max :      | 14.93 | kN    | Punch NO.; | 45 |
| Filling Depth :      | Increase ; 0 | Times | Decrease ; | 15    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.10    | mm  | Main Thickness ; | 2.11     | mm           |
| Filling Depth :               | 6.67    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.742   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 14 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 21 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 550304             | Tablets |       |        |         |
| A.W.C Counter : | 545613             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 545613<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | vg96371       | Batch NO. :          | VIM2500340            |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 33 Min 56 Sec |
| Total Counter :  | 550304        | Total Running Time : | 514 Hour 0 Min        |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |            |
|------------------------------------|-------------|----------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 76       | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm    |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |            |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |            |
| Oil Lubrication :                  |             |          |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec        |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |            |
|                                    |             |          |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 20.07    | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 19.30    | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 15.75    | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 15.44    | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 15.13    | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 11.58    | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 10.81    | kN      | No                 |              |            |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |            |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |          |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.16       | kN | Min : | 1.99       | kN | Punch NO.; | 7  |
| SD :                 | 3.20       | %  | Max : | 2.33       | kN | Punch NO.; | 46 |
| Mean Main Pressure : | 15.29      | kN | Min : | 14.78      | kN | Punch NO.; | 33 |
| SD :                 | 1.50       | %  | Max : | 15.83      | kN | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 0  | Times | Decrease ; | 1  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.20    | mm  | Main Thickness ; | 2.21     | mm           |
| Filling Depth :               | 7.62    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.742   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 6  |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 21 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 550304             | Tablets |       |        |         |
| A.W.C Counter : | 545613             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 545613<br>(100.0%) | Tablets |       |        |         |