

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | vg96371       | Batch NO. :          | VIM2500340            |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 12 Min 32 Sec |
| Total Counter :  | 650352        | Total Running Time : | 514 Hour 38 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 75        | %       | Manual ;           | 35    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 1500000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 18.87     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 18.14     | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 14.80     | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 14.51     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 14.22     | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.88     | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 10.16     | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |           |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.13       | kN | Min : | 2.04       | kN | Punch NO.; | 24 |
| SD :                 | 2.20       | %  | Max : | 2.25       | kN | Punch NO.; | 6  |
| Mean Main Pressure : | 14.46      | kN | Min : | 14.13      | kN | Punch NO.; | 15 |
| SD :                 | 1.20       | %  | Max : | 14.92      | kN | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 0  | Times | Decrease ; | 15 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.10    | mm  | Main Thickness ; | 2.11     | mm           |
| Filling Depth :               | 6.67    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 15      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 14 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 21 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 650352             | Tablets |       |        |         |
| A.W.C Counter : | 645648             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 645648<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | ACICLO400MG,STV11B0   |
| User ID :        | vg96371       | Batch NO. :          | VIM2500340            |
| Print Interval : | 50000 Tablets | Running Time :       | 04 Hour 12 Min 32 Sec |
| Total Counter :  | 650352        | Total Running Time : | 514 Hour 38 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 76        | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 1500000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 20.07     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 19.30     | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 15.75     | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 15.44     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 15.13     | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 11.58     | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 10.81     | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.12       | kN | Min : | 1.98       | kN | Punch NO.; | 41 |
| SD :                 | 3.20       | %  | Max : | 2.27       | kN | Punch NO.; | 51 |
| Mean Main Pressure : | 15.21      | kN | Min : | 14.52      | kN | Punch NO.; | 33 |
| SD :                 | 1.70       | %  | Max : | 15.70      | kN | Punch NO.; | 18 |
| Filling Depth :      | Increase ; | 0  | Times | Decrease ; | 1  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 156000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.20    | mm  | Main Thickness ; | 2.21     | mm           |
| Filling Depth :               | 7.62    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 15      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.01 | kN | Max Punch ; | 14 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 14 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 21 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 650352             | Tablets |       |        |         |
| A.W.C Counter : | 645648             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 645648<br>(100.0%) | Tablets |       |        |         |