

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pm98770	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 39 Sec
Total Counter :	51392	Total Running Time :	576 Hour 31 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	150	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		600000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	23.13	kN	No		
HEP :	25.0 %	22.24	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	18.15	kN	Yes 15	Times	
Ref :		17.79	kN			
LCP :	2.0 %	17.43	kN	Yes 15	Times	
LEP :	25.0 %	13.34	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	12.45	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	2.88	kN	Min :	2.47	kN	Punch NO.;	13
		7.90	%	Max :	3.61	kN	Punch NO.;	42
Mean Main Pressure :		18.39	kN	Min :	16.19	kN	Punch NO.;	13
	SD :	5.40	%	Max :	21.96	kN	Punch NO.;	42
Filling Depth :	Increase ;	20	Times	Decrease ;	5	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.04	mm	Main Thickness ;	1.05	mm
Filling Depth :	8.27	mm	Current Cam :	Cam C	
Main Air pressure :	0.809	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	1	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.39	kN		
	Max ;	0.49	kN	Max Punch ;	52
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.30	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	51392	Tablets			
A.W.C Counter :	44636	Tablets			
HEP :	37 (0.1%)	Tablets	LEP :	32 (0.1%)	Tablets
GOOD :	44567 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pm98770	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	00 Hour 32 Min 39 Sec
Total Counter :	51392	Total Running Time :	576 Hour 31 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	150	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		600000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	25.06	kN	No		
HEP :	25.0 %	24.10	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	19.67	kN	Yes 15	Times	
Ref :		19.28	kN			
LCP :	2.0 %	18.89	kN	Yes 15	Times	
LEP :	25.0 %	14.46	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	13.50	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.79	kN	Min :	2.41	kN	Punch NO.;	40
		7.10	%	Max :	3.38	kN	Punch NO.;	30
Mean Main Pressure :		19.65	kN	Min :	17.83	kN	Punch NO.;	40
	SD :	4.80	%	Max :	21.95	kN	Punch NO.;	30
Filling Depth :	Increase ;	15	Times	Decrease ;	5	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	4.12	mm	Main Thickness ;	1.12	mm
Filling Depth :	8.02	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.809	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	1	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.39	kN		
	Max ;	0.49	kN	Max Punch ;	52
Lower Punch :	Average ;	0.19	kN		
	Max ;	0.30	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	51392	Tablets			
A.W.C Counter :	44636	Tablets			
HEP :	37 (0.1%)	Tablets	LEP :	32 (0.1%)	Tablets
GOOD :	44567 (99.8%)	Tablets			