

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pb97604	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	00 Hour 56 Min 38 Sec
Total Counter :	101388	Total Running Time :	576 Hour 55 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	150	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		600000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	22.91	kN	No		
HEP :	25.0 %	22.03	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	17.98	kN	Yes 15	Times	
Ref :		17.62	kN			
LCP :	2.0 %	17.27	kN	Yes 15	Times	
LEP :	25.0 %	13.22	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	12.34	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.71	kN	Min :	2.29	kN	Punch NO.;	27
SD :	9.00	%	Max :	3.38	kN	Punch NO.;	46
Mean Main Pressure :	17.60	kN	Min :	15.21	kN	Punch NO.;	3
SD :	7.10	%	Max :	20.39	kN	Punch NO.;	46
Filling Depth :	Increase ; 36	Times	Decrease ;	20	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	30	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.04	mm	Main Thickness ;	1.05	mm
Filling Depth :	8.55	mm	Current Cam :	Cam C	
Main Air pressure :	0.809	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.14	kN		
	Max ;	0.24	kN	Max Punch ;	50
Lower Punch :	Average ;	0.20	kN		
	Max ;	0.32	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	101388	Tablets			
A.W.C Counter :	88750	Tablets			
HEP :	44 (0.0%)	Tablets	LEP :	32 (0.0%)	Tablets
GOOD :	88674 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	Pentoxifylline 400mg
User ID :	pb97604	Batch NO. :	PEN TAB XV 25 001
Print Interval :	50000 Tablets	Running Time :	00 Hour 56 Min 38 Sec
Total Counter :	101388	Total Running Time :	576 Hour 55 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	150	%	Manual ;	15	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		600000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	10.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	23.98	kN	No		
HEP :	25.0 %	23.06	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	18.81	kN	Yes	15	Times
Ref :		18.44	kN			
LCP :	2.0 %	18.08	kN	Yes	15	Times
LEP :	25.0 %	13.83	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	12.91	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.59	kN	Min :	2.09	kN	Punch NO.;	23
		7.20	%	Max :	3.05	kN	Punch NO.;	30
Mean Main Pressure :		18.76	kN	Min :	16.39	kN	Punch NO.;	23
	SD :	4.90	%	Max :	20.76	kN	Punch NO.;	30
Filling Depth :	Increase ;	26	Times	Decrease ;	19	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	30	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	4.80	mm	Pre Thickness ;	1.81	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	4.12	mm	Main Thickness ;	1.12	mm
Filling Depth :	8.02	mm	Current Cam :	Cam C	
Main Air pressure :	0.809	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	11	Min			
Lower Punch Remain Time (S2):	14	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.14	kN		
	Max ;	0.24	kN	Max Punch ;	50
Lower Punch :	Average ;	0.20	kN		
	Max ;	0.33	kN	Max Punch ;	15
Ejection Force :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0

Tablet Data (N/A)

Total Counter :	101388	Tablets			
A.W.C Counter :	88750	Tablets			
HEP :	44 (0.0%)	Tablets	LEP :	32 (0.0%)	Tablets
GOOD :	88674 (99.9%)	Tablets			