

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | yu98249       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 21 Hour 53 Min 16 Sec   |
| Total Counter :  | 5668848       | Total Running Time : | 631 Hour 51 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 22                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 18.08             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 17.38             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 14.18             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 13.90             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 13.63             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 10.43             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 9.73              | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.08       | kN  | Min : | 1.03       | kN  | Punch NO.; | 22 |
| SD :                 | 2.30       | %   | Max : | 1.17       | kN  | Punch NO.; | 75 |
| Mean Main Pressure : | 14.03      | kN  | Min : | 12.89      | kN  | Punch NO.; | 22 |
| SD :                 | 4.30       | %   | Max : | 15.98      | kN  | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 386 | Times | Decrease ; | 423 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.80   | mm  | Pre Thickness ;  | 2.81     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.51   | mm  | Main Thickness ; | 1.52     | mm           |
| Filling Depth :               | 5.20   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.800  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 27     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.28 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 35 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 5668848          | Tablets |       |           |         |
| A.W.C Counter : | 5664267          | Tablets |       |           |         |
| HEP :           | 204 (0.0%)       | Tablets | LEP : | 99 (0.0%) | Tablets |
| GOOD :          | 5663964 (100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SETRA 50<br>MG,STSH4A00 |
| User ID :        | yu98249       | Batch NO. :          | SHM2500350              |
| Print Interval : | 50000 Tablets | Running Time :       | 21 Hour 53 Min 16 Sec   |
| Total Counter :  | 5668848       | Total Running Time : | 631 Hour 51 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 25                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 10.0  | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 1.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 19.10             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 18.36             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 14.99             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 14.69             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 14.40             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 11.02             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 10.28             | kN      | No                 |       |                      |
| SD limit :                         | 20.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 20.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.20       | kN  | Min : | 1.15       | kN  | Punch NO.; | 77 |
| SD :                 | 2.30       | %   | Max : | 1.33       | kN  | Punch NO.; | 35 |
| Mean Main Pressure : | 14.67      | kN  | Min : | 13.39      | kN  | Punch NO.; | 5  |
| SD :                 | 4.80       | %   | Max : | 17.61      | kN  | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 182 | Times | Decrease ; | 220 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.70    | mm  | Pre Thickness ;  | 2.71     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.64    | mm  | Main Thickness ; | 1.64     | mm           |
| Filling Depth :               | 5.63    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.800   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 27      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 35 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 5668848             | Tablets |       |           |         |
| A.W.C Counter : | 5664267             | Tablets |       |           |         |
| HEP :           | 204 (0.0%)          | Tablets | LEP : | 99 (0.0%) | Tablets |
| GOOD :          | 5663964<br>(100.0%) | Tablets |       |           |         |