

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET1000mg STYQ1A00    |
| User ID :        | mr99799       | Batch NO. :          | YQM2501110            |
| Print Interval : | 50000 Tablets | Running Time :       | 16 Hour 57 Min 35 Sec |
| Total Counter :  | 2150916       | Total Running Time : | 691 Hour 8 Min        |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 120      | %       | Manual ;           | 40           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 27.63    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 26.57    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 21.68    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 21.25    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 20.83    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 15.94    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 14.88    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.13  | kN    | Min :      | 2.86  | kN    | Punch NO.; | 3  |
|                      |            | 3.10  | %     | Max :      | 3.35  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : |            | 21.13 | kN    | Min :      | 19.51 | kN    | Punch NO.; | 3  |
|                      | SD :       | 2.60  | %     | Max :      | 22.47 | kN    | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 86    | Times | Decrease ; | 115   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.70   | mm  | Pre Thickness ;  | 4.71     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.64   | mm  | Main Thickness ; | 3.65     | mm           |
| Filling Depth :               | 9.19   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.799  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 45 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 7  |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2150916             | Tablets |       |           |         |
| A.W.C Counter : | 2148484             | Tablets |       |           |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 64 (0.0%) | Tablets |
| GOOD :          | 2148418<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET1000mg STYQ1A00    |
| User ID :        | mr99799       | Batch NO. :          | YQM2501110            |
| Print Interval : | 50000 Tablets | Running Time :       | 16 Hour 57 Min 35 Sec |
| Total Counter :  | 2150916       | Total Running Time : | 691 Hour 8 Min        |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 120      | %       | Manual ;           | 40           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.67    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 28.53    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 23.28    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 22.83    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 22.37    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 17.12    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 15.98    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 3.12          | kN    | Min :         | 2.95  | kN | Punch NO.; | 3  |
| SD :                 | 2.40          | %     | Max :         | 3.30  | kN | Punch NO.; | 8  |
| Mean Main Pressure : | 22.91         | kN    | Min :         | 22.01 | kN | Punch NO.; | 28 |
| SD :                 | 1.90          | %     | Max :         | 23.79 | kN | Punch NO.; | 17 |
| Filling Depth :      | Increase ; 49 | Times | Decrease ; 58 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20     | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.70   | mm  | Pre Thickness ;  | 4.71     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.64   | mm  | Main Thickness ; | 3.65     | mm           |
| Filling Depth :               | 9.40   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.799  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 45 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 7  |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2150916             | Tablets |       |           |         |
| A.W.C Counter : | 2148484             | Tablets |       |           |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 64 (0.0%) | Tablets |
| GOOD :          | 2148418<br>(100.0%) | Tablets |       |           |         |