

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	yu98249	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	00 Hour 18 Min 06 Sec
Total Counter :	28502	Total Running Time :	710 Hour 12 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	141	%	Manual ;	66	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2025000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	20.12	kN	No		
HEP :	25.0 %	19.35	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	15.79	kN	Yes 15	Times	
Ref :		15.48	kN			
LCP :	2.0 %	15.17	kN	Yes 15	Times	
LEP :	25.0 %	11.61	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	10.84	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	0.47	kN	Min :	0.44	kN	Punch NO.;	9
SD :	3.10	%	Max :	0.52	kN	Punch NO.;	42
Mean Main Pressure :	15.06	kN	Min :	13.37	kN	Punch NO.;	15
SD :	4.70	%	Max :	17.70	kN	Punch NO.;	42
Filling Depth :	Increase ;	329	Times	Decrease ;	320	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	21	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	7.00	mm	Pre Thickness ;	4.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.62	mm	Main Thickness ;	1.63	mm
Filling Depth :	6.81	mm	Current Cam :	Cam B	
Main Air pressure :	0.799	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	47	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.12	kN	Max Punch ;	43
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.09	kN	Max Punch ;	4
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	2

Tablet Data (N/A)

Total Counter :	28502	Tablets			
A.W.C Counter :	25138	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	25138 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	yu98249	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	00 Hour 18 Min 06 Sec
Total Counter :	28502	Total Running Time :	710 Hour 12 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	149	%	Manual ;	74	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		2025000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	20.33	kN	No		
HEP :	25.0 %	19.55	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	15.95	kN	Yes 15	Times	
Ref :		15.64	kN			
LCP :	2.0 %	15.33	kN	Yes 15	Times	
LEP :	25.0 %	11.73	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	10.95	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	0.64	kN	Min :	0.60	kN	Punch NO.;	24
SD :	1.90	%	Max :	0.67	kN	Punch NO.;	2
Mean Main Pressure :	15.51	kN	Min :	14.76	kN	Punch NO.;	21
SD :	2.60	%	Max :	16.62	kN	Punch NO.;	20
Filling Depth :	Increase ;	157	Times	Decrease ;	157	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	22	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	7.00	mm	Pre Thickness ;	4.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.70	mm	Main Thickness ;	1.71	mm
Filling Depth :	6.49	mm	Current Cam :	Cam B	
Main Air pressure :	0.799	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	47	Min			
Lower Punch Remain Time (S2):	17	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.06	kN		
	Max ;	0.12	kN	Max Punch ;	43
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.09	kN	Max Punch ;	4
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.05	kN	Max Punch ;	2

Tablet Data (N/A)

Total Counter :	28502	Tablets			
A.W.C Counter :	25138	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	25138 (100.0%)	Tablets			