

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	kp97154	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	10 Hour 01 Min 49 Sec
Total Counter :	1450548	Total Running Time :	719 Hour 56 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	105	%	Manual ;	66	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		2.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	25.69	kN	No		
HEP :	25.0 %	24.70	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	20.15	kN	Yes 15	Times	
Ref :		19.76	kN			
LCP :	2.0 %	19.36	kN	Yes 15	Times	
LEP :	25.0 %	14.82	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	13.83	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	SD :	0.63	kN	Min :	0.57	kN	Punch NO.;	21
		4.00	%	Max :	0.72	kN	Punch NO.;	19
Mean Main Pressure :	SD :	19.92	kN	Min :	18.09	kN	Punch NO.;	21
		4.70	%	Max :	22.23	kN	Punch NO.;	15
Filling Depth :	Increase ;	1260	Times	Decrease ;	1274	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	21	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.70	mm	Pre Thickness ;	3.71	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.50	mm	Main Thickness ;	1.51	mm
Filling Depth :	6.84	mm	Current Cam :	Cam C	
Main Air pressure :	0.722	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	28	Min			
Lower Punch Remain Time (S2):	28	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.21	kN		
	Max ;	1.80	kN	Max Punch ;	18
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	1450548	Tablets			
A.W.C Counter :	1427152	Tablets			
HEP :	1039 (0.1%)	Tablets	LEP :	651 (0.0%)	Tablets
GOOD :	1425462 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO400MG STIU2B00
User ID :	kp97154	Batch NO. :	IUM2500140
Print Interval :	50000 Tablets	Running Time :	10 Hour 01 Min 49 Sec
Total Counter :	1450548	Total Running Time :	719 Hour 56 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	126	%	Manual ;	74	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		30	Min	Supply Time ;	0.5	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		2.00	kN			
Ejecting Force High Limit :		1.00	kN			
				Stop Condition ;		
HSP :	30.0 %	27.12	kN	No		
HEP :	25.0 %	26.08	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	21.28	kN	Yes 15	Times	
Ref :		20.86	kN			
LCP :	2.0 %	20.44	kN	Yes 15	Times	
LEP :	25.0 %	15.65	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	14.60	kN	No		
SD limit :	20.0 %			Yes		
Pre HSP		20.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	0.79	kN	Min :	0.74	kN	Punch NO.;	14
SD :	3.40	%	Max :	0.94	kN	Punch NO.;	19
Mean Main Pressure :	20.66	kN	Min :	18.59	kN	Punch NO.;	18
SD :	3.10	%	Max :	22.13	kN	Punch NO.;	33
Filling Depth :	Increase ;	812	Times	Decrease ;	772	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	25	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.70	mm	Pre Thickness ;	3.71	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.55	mm	Main Thickness ;	1.56	mm
Filling Depth :	7.22	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.722	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	28	Min			
Lower Punch Remain Time (S2):	28	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.21	kN		
	Max ;	1.80	kN	Max Punch ;	18
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	1450548	Tablets			
A.W.C Counter :	1427205	Tablets			
HEP :	1039 (0.1%)	Tablets	LEP :	651 (0.0%)	Tablets
GOOD :	1425515 (99.9%)	Tablets			