

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAl3A00  |
| User ID :        | am98755       | Batch NO. :          | AIM2500220            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 08 Min 37 Sec |
| Total Counter :  | 800630        | Total Running Time : | 731 Hour 20 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |       |                      |
|------------------------------------|-------------|----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 105      | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min : | 2.50 mm              |
|                                    |             | 4.5-10mm |         |                    | Max : | 8.00 mm              |
| Target Quantity :                  |             | 1500000  | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |          |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5   | Sec                  |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 2        | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |       |                      |
|                                    |             |          |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 25.60    | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 24.62    | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.5 %       | 20.19    | kN      | Yes                | 15    | Times                |
| Ref :                              |             | 19.70    | kN      |                    |       |                      |
| LCP :                              | 2.5 %       | 19.20    | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 %      | 14.77    | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 %      | 13.79    | kN      | No                 |       |                      |
| SD limit :                         | 20.0 %      |          |         | Yes                |       |                      |
| Pre HSP                            |             | 20.00    | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None        |          |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.69       | kN  | Min : | 2.49       | kN  | Punch NO.; | 10 |
| SD :                 | 4.00       | %   | Max : | 2.92       | kN  | Punch NO.; | 46 |
| Mean Main Pressure : | 18.63      | kN  | Min : | 17.48      | kN  | Punch NO.; | 16 |
| SD :                 | 3.20       | %   | Max : | 20.04      | kN  | Punch NO.; | 46 |
| Filling Depth :      | Increase ; | 280 | Times | Decrease ; | 299 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 17     | RPM | Capacity :       | 108120   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.69   | mm  | Pre Thickness ;  | 1.70     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.83   | mm  | Main Thickness ; | 0.84     | mm           |
| Filling Depth :               | 5.70   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.788  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 11     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | On       |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 800630            | Tablets |       |            |         |
| A.W.C Counter : | 783130            | Tablets |       |            |         |
| HEP :           | 715 (0.1%)        | Tablets | LEP : | 558 (0.1%) | Tablets |
| GOOD :          | 781857<br>(99.8%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amioda 200mgSTAl3A00  |
| User ID :        | am98755       | Batch NO. :          | AIM2500220            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 08 Min 37 Sec |
| Total Counter :  | 800630        | Total Running Time : | 731 Hour 20 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 96       | %       | Manual ;           | 32           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 0.5          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 1.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 25.47    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 24.49    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.5 %       | 20.08    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 19.59    | kN      |                    |              |           |
| LCP :                              | 2.5 %       | 19.10    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 14.70    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 13.72    | kN      | No                 |              |           |
| SD limit :                         | 20.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 20.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.78  | kN    | Min :      | 2.49  | kN    | Punch NO.; | 24 |
|                      |            | 5.30  | %     | Max :      | 3.18  | kN    | Punch NO.; | 17 |
| Mean Main Pressure : |            | 19.35 | kN    | Min :      | 18.15 | kN    | Punch NO.; | 24 |
|                      | SD :       | 3.30  | %     | Max :      | 21.03 | kN    | Punch NO.; | 45 |
| Filling Depth :      | Increase ; | 271   | Times | Decrease ; | 260   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 17      | RPM | Capacity :       | 108120   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.70    | mm  | Pre Thickness ;  | 1.71     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.90    | mm  | Main Thickness ; | 0.91     | mm           |
| Filling Depth :               | 6.26    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.788   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 1       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 11      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | On       |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 17 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 800630            | Tablets |       |            |         |
| A.W.C Counter : | 783130            | Tablets |       |            |         |
| HEP :           | 715 (0.1%)        | Tablets | LEP : | 558 (0.1%) | Tablets |
| GOOD :          | 781857<br>(99.8%) | Tablets |       |            |         |