

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | sj97904       | Batch NO. :          | SHM2500910            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 58 Min 23 Sec |
| Total Counter :  | 250230        | Total Running Time : | 801 Hour 14 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 26       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60       | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 18.33    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 17.62    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 14.38    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 14.10    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 13.82    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 10.57    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 9.87     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.56  | kN    | Min :      | 2.44  | kN    | Punch NO.; | 16 |
|                      |            | 2.20  | %     | Max :      | 2.71  | kN    | Punch NO.; | 71 |
| Mean Main Pressure : |            | 14.05 | kN    | Min :      | 13.28 | kN    | Punch NO.; | 43 |
|                      | SD :       | 2.60  | %     | Max :      | 14.81 | kN    | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 20    | Times | Decrease ; | 21    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 3.00     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.37   | mm  | Main Thickness ; | 1.87     | mm           |
| Filling Depth :               | 7.76   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.788  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.38 | kN |             |    |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 23 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 58 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 74 |

### Tablet Data (N/A)

|                 |          |         |       |        |         |
|-----------------|----------|---------|-------|--------|---------|
| Total Counter : | 250230   | Tablets |       |        |         |
| A.W.C Counter : | 243883   | Tablets |       |        |         |
| HEP :           | 1 (0.0%) | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 243882   | Tablets |       |        |         |
|                 | (100.0%) |         |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | sj97904       | Batch NO. :          | SHM2500910            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 58 Min 23 Sec |
| Total Counter :  | 250230        | Total Running Time : | 801 Hour 14 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 25       | %       | Manual ;           | 48           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60       | Min     | Supply Time ;      | 0.3          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 20.31    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 19.53    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 15.94    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 15.63    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 15.31    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 11.72    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 10.94    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.62          | kN    | Min :      | 2.51  | kN    | Punch NO.; | 49 |
| SD :                 | 1.90          | %     | Max :      | 2.73  | kN    | Punch NO.; | 42 |
| Mean Main Pressure : | 15.54         | kN    | Min :      | 14.75 | kN    | Punch NO.; | 11 |
| SD :                 | 2.40          | %     | Max :      | 16.41 | kN    | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 22 | Times | Decrease ; | 38    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 3.00     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.43   | mm  | Main Thickness ; | 1.93     | mm           |
| Filling Depth :               | 7.48   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.788  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.38 | kN |             |    |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 23 |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 58 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 74 |

### Tablet Data (N/A)

|                 |          |         |       |        |         |
|-----------------|----------|---------|-------|--------|---------|
| Total Counter : | 250230   | Tablets |       |        |         |
| A.W.C Counter : | 243883   | Tablets |       |        |         |
| HEP :           | 1 (0.0%) | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 243882   | Tablets |       |        |         |
|                 | (100.0%) |         |       |        |         |