

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	pb97604	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	05 Hour 29 Min 55 Sec
Total Counter :	694872	Total Running Time :	817 Hour 17 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	80	%	Manual ;	25	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.62	kN	No		
HEP :	25.0 %	16.94	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	13.83	kN	Yes 10	Times	
Ref :		13.55	kN			
LCP :	2.0 %	13.28	kN	Yes 10	Times	
LEP :	25.0 %	10.17	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	9.49	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.13	kN	Min :	2.07	kN	Punch NO.;	28
SD :	2.50	%	Max :	2.46	kN	Punch NO.;	48
Mean Main Pressure :	13.72	kN	Min :	12.46	kN	Punch NO.;	19
SD :	5.00	%	Max :	15.66	kN	Punch NO.;	3
Filling Depth :	Increase ;	215	Times	Decrease ;	177	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	16	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.65	mm	Main Thickness ;	1.66	mm
Filling Depth :	9.48	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.820	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	4	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.16	kN		
	Max ;	0.25	kN	Max Punch ;	42
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.52	kN	Max Punch ;	48
Ejection Force :	Average ;	0.10	kN		
	Max ;	1.02	kN	Max Punch ;	48

Tablet Data (N/A)

Total Counter :	694872	Tablets			
A.W.C Counter :	689660	Tablets			
HEP :	123 (0.0%)	Tablets	LEP :	18 (0.0%)	Tablets
GOOD :	689519 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	pb97604	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	05 Hour 29 Min 55 Sec
Total Counter :	694872	Total Running Time :	817 Hour 17 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	80	%	Manual ;	25	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.45	kN	No		
HEP :	25.0 %	16.78	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	13.69	kN	Yes 10	Times	
Ref :		13.43	kN			
LCP :	2.0 %	13.16	kN	Yes 10	Times	
LEP :	25.0 %	10.07	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	9.40	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.04	kN	Min :	1.98	kN	Punch NO.;	28
SD :	2.30	%	Max :	2.33	kN	Punch NO.;	48
Mean Main Pressure :	13.35	kN	Min :	12.22	kN	Punch NO.;	48
SD :	4.40	%	Max :	14.55	kN	Punch NO.;	43
Filling Depth :	Increase ;	51	Times	Decrease ;	40	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	16	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.65	mm	Main Thickness ;	1.66	mm
Filling Depth :	9.27	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.820	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	4	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.16	kN		
	Max ;	0.24	kN	Max Punch ;	41
Lower Punch :	Average ;	0.11	kN		
	Max ;	0.52	kN	Max Punch ;	48
Ejection Force :	Average ;	0.10	kN		
	Max ;	1.02	kN	Max Punch ;	48

Tablet Data (N/A)

Total Counter :	694872	Tablets			
A.W.C Counter :	689660	Tablets			
HEP :	123 (0.0%)	Tablets	LEP :	18 (0.0%)	Tablets
GOOD :	689519 (100.0%)	Tablets			