

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	sj97904	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	14 Hour 57 Min 33 Sec
Total Counter :	1895244	Total Running Time :	826 Hour 46 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	102	%	Manual ;	25	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2000000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	18.05	kN	No		
HEP :	25.0 %	17.35	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	14.16	kN	Yes 10	Times	
Ref :		13.88	kN			
LCP :	2.0 %	13.61	kN	Yes 10	Times	
LEP :	25.0 %	10.41	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	9.72	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.31	kN	Min :	2.23	kN	Punch NO.;	49
SD :	3.30	%	Max :	2.63	kN	Punch NO.;	8
Mean Main Pressure :	13.55	kN	Min :	12.44	kN	Punch NO.;	26
SD :	3.80	%	Max :	14.60	kN	Punch NO.;	8
Filling Depth :	Increase ;	559	Times	Decrease ;	557	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	20	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.68	mm	Main Thickness ;	1.69	mm
Filling Depth :	9.02	mm	Current Cam :	Cam C	
Main Air pressure :	0.791	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	18	Min			
Lower Punch Remain Time (S2):	19	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.21	kN		
	Max ;	0.28	kN	Max Punch ;	41
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.56	kN	Max Punch ;	18
Ejection Force :	Average ;	0.15	kN		
	Max ;	0.94	kN	Max Punch ;	18

Tablet Data (N/A)

Total Counter :	1895244	Tablets			
A.W.C Counter :	1889126	Tablets			
HEP :	636 (0.0%)	Tablets	LEP :	326 (0.0%)	Tablets
GOOD :	1888164 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MEBEV 135 mgSTMV2A00
User ID :	sj97904	Batch NO. :	MVM2500413
Print Interval :	50000 Tablets	Running Time :	14 Hour 57 Min 33 Sec
Total Counter :	1895244	Total Running Time :	826 Hour 46 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	102	%	Manual ;	25	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2000000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		30	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		20	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.38	kN	No		
HEP :	25.0 %	16.71	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	13.64	kN	Yes 10	Times	
Ref :		13.37	kN			
LCP :	2.0 %	13.10	kN	Yes 10	Times	
LEP :	25.0 %	10.03	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	9.36	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.35	kN	Min :	2.26	kN	Punch NO.;	18
SD :	1.60	%	Max :	2.50	kN	Punch NO.;	36
Mean Main Pressure :	13.24	kN	Min :	11.51	kN	Punch NO.;	19
SD :	4.20	%	Max :	14.36	kN	Punch NO.;	15
Filling Depth :	Increase ;	454	Times	Decrease ;	464	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	127200	Tablets/hour
Feeder : Status ;	Auto		Speed ;	20	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.69	mm	Main Thickness ;	1.70	mm
Filling Depth :	8.96	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.791	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	18	Min			
Lower Punch Remain Time (S2):	19	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.21	kN		
	Max ;	0.28	kN	Max Punch ;	41
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.56	kN	Max Punch ;	18
Ejection Force :	Average ;	0.15	kN		
	Max ;	0.92	kN	Max Punch ;	18

Tablet Data (N/A)

Total Counter :	1895244	Tablets			
A.W.C Counter :	1889126	Tablets			
HEP :	636 (0.0%)	Tablets	LEP :	326 (0.0%)	Tablets
GOOD :	1888164 (99.9%)	Tablets			